



YAFHO BIO-TECHNOLOGY CO., LTD.


Tel: 86-20-85295649 Fax: 86-20-85295797
<http://www.yafho.com> E-mail: info@yafho.com

YAFHO BIO-TECHNOLOGY CO., LTD.



Yafho is one of the professional wound dressing **Manufacturers** with more than **10000 m2 GMP standard workshop**.

With scientific and rigorous management and excellent team, we can ensure the product quality. We sincerely look forward to cooperate with you.



Rev.:YE-E201909

Sili- Care

Sili-Care, which consists of a gentle silicone wound contact layer, a super absorbent pad layer, and a vapor permeable and waterproof film, is an ideal dressing for advanced wound exudate management. The multi-layers construction facilitates dynamic fluid management to provide optimal moist wound environment, which promotes wound healing, minimizes trauma and reduces pain during dressing changes.


Benefits

- Conformable dressing – provide super gentle silicone contact layer for fragile skin.
- Minimize trauma and pain during dressing change.
- Easy to apply and remove – the dressing can be repositioned easily without losing its adherent properties.
- Shower proof dressing allowing patients to shower with the dressing in place.
- Customized shape design to adapt different requirements.
- Maintains its integrity during use

Usage:

- Pressure ulcer
- Leg ulcer
- Non-infected diabetic foot ulcer
- Donor sites
- Skin abrasions
- Post-operative wound dressing
- Lacerations and abrasion

Sili cone
Dressing



SILI-CARE (Border)	SILI CARE (Border Lite)	Size(cm)	Package
626375	626475	7.5x7.5(4.5x4.5)	10pcs/box,12boxes/ctn
626310	626410	10x10(6x6)	10pcs/box,12boxes/ctn
626313	626413	12.5x12.5(8.5x8.5)	10pcs/box,6boxes/ctn
626315	626415	10x15(6x11)	10pcs/box,6boxes/ctn
626355	626455	15x15(11x11)	10pcs/box,6boxes/ctn
626320	626420	10x20(6x16)	10pcs/box,6boxes/ctn
626352	626452	15x20(11x16)	10pcs/box,6boxes/ctn
626325	626425	10x25(6x21)	10pcs/box,6boxes/ctn
626330	626430	10x30(6x26)	10pcs/box,6boxes/ctn
626318	626418	18x18	10pcs /box,6boxes /ctn
626323	626423	23x23	10pcs /box,6boxes /ctn

SILI-CARE	SILI-CARE Lite	Size(cm)	Package
626210	626110	10x10	10pcs /box,12boxes /ctn
626213	626113	12.5x12.5	10pcs /box,6boxes /ctn
626255	626155	15x15	10pcs /box,6boxes /ctn
626220	626120	10x20	10pcs /box,6boxes /ctn
626225	626125	10x25	10pcs /box,6boxes /ctn
626230	626130	10x30	10pcs /box,6boxes /ctn

Tel:86-20-85295649 Fax:86-20-85295797 <http://www.yafho.com> E-mail:info@yafho.com

Foam Dressing

Foam dressing is a highly absorbent dressing made from hydrophilic polyurethane foam, ideally designed for heavily exuding wounds. The soft flexible foam absorbs large amount drainage, ensures patients comfort; the top backing film protect wounds against external bacterial.

Applications:


Foam dressing is indicated for the moderate to heavily exuding wounds
Pressure ulcer
Leg venous ulcer
Diabetic foot ulcer
First and second degree burns
Skin grafts and donor site
Post-operative and traumatic wounds





Ref.	Size(cm)	Package
628001	7.5x7.5	10pcs/box,20boxes/ctn
628002	10x10	10pcs/box,20boxes/ctn
628003	12.5x12.5	10pcs/box,10boxes/ctn
628004	15x15	10pcs/box,10boxes/ctn
628005	15x20	10pcs/box,10boxes/ctn

Silicone Non-woven Dressing

Silicone nonwoven wound dressing is updated with silicone glue. It can efficiently protect against tape-related skin injury. It is specially designed for patients injuries require debridement, wound irrigation, mobility treatment, cushioning, antimicrobials or surgical grafting.




Ref.	Size(cm)	Package
616105	5x8(2.2x4.5)	100pcsx20boxes,2000pcs/ctn
616106	6x7(3x4)	100pcsx20boxes,2000pcs/ctn
616108	6x10(3x7)	50pcsx24boxes,1200pcs/ctn
616110	9x10(5x6)	30pcsx24boxes,720pcs/ctn
616115	9x15(5x10)	30pcsx24boxes,720pcs/ctn
616120	10x20(5x15)	25pcsx18boxes,450pcs/ctn
616125	10x25(5x20)	25pcsx18boxes,450pcs/ctn



Yafho
WOUND CARE

Negative Pressure Wound Therapy System




NPWT Single Use Negative Pressure Wound Therapy system is indicated for patients who would benefit from a suction device (Negative Pressure Wound Therapy) as it may promote wound healing via removal of low to moderate levels of exudate and infectious materials. NPWT Single Use Negative Pressure Wound Therapy System is suitable for use in both a hospital and community setting and approved for the following indications.

- Acute wounds
- Chronic wounds
- Flaps and grafts
- Partial-thickness burns
- Subacute and dehiscent wounds
- Surgically closed incision sites
- Traumatic wounds

Vacuum Sealing Drainage(VSD)Dressing Set



The VSD dressings are designed to heal Hard-to-Heal wounds, such as Diabetic Foot, Decubitus, Burn Wounds, etc. Composition: Medical Sponge, Film Dressings, Drainage Tubes, Clamps.



Yafho
WOUND CARE


Wound Dressing

Island Non woven Dressing

The dressing is made of non-woven fabric and non-adherent pad. It is used for wound care.

Features	Normal	Superior
Water-repellent		★★★★★
Breathable	★★★	★★★★★
Adhesive	★★★	★★★★★
Hypoallergenic	★★★★★	★★★★★
Absorbent	★★★★★	★★★★★

Ref. (Normal)	Ref. (Superior)	Size(cm)	Package
611105	610105	5x8(2.2x4.5)	100pcsx20boxes,2000pcs/ctn
611106	610106	6x7(3x4)	100pcsx20boxes,2000pcs/ctn
611108	610108	6x10(3x7)	50pcsx24boxes,1200pcs/ctn
611110	610110	9x10(5x6)	30pcsx24boxes,720pcs/ctn
611115	610115	9x15(5x10)	30pcsx24boxes,720pcs/ctn
611120	610120	10x20(5x15)	25pcsx18boxes,450pcs/ctn
611125	610125	10x25(5x20)	25pcsx18boxes,450pcs/ctn
611130	610130	10x30(5x25)	25pcsx12boxes,300pcs/ctn
611135	610135	10x35(5x30)	25pcsx12boxes,300pcs/ctn
611140	610140	10x40(5x35)	25pcsx12boxes,300pcs/ctn
611065	610065	8x8(φ 5)	100pcsx18boxes,1800pcs/ctn
611086	610086	12x12(φ 8)	100pcsx18boxes,1800pcs/ctn
611087	610087	7.5x7.5(φ 4.5)	100pcsx18boxes,1800pcs/ctn
611024	610124	12x14(8x9.5)	30pcsx24boxes,720pcs/ctn
611012	610112	9x10(5x6)	50pcsx24boxes,1200pcs/ctn



Non – adherent Dressing Pad

Ref.	Size(cm)	Package
688001	5x5	50pcsx36boxes,1800pcs/ctn
688002	7.5x7.5	50pcsx36boxes,1800pcs/ctn
688003	5x7.5	50pcsx36boxes,1800pcs/ctn
688004	10x7.5	50pcsx27boxes,1350pcs/ctn
688005	10x10	50pcsx27boxes,1350pcs/ctn



Yafho
WOUND CARE

Tel:86-20-85295649 Fax:86-20-85295797 <http://www.yafho.com> E-mail:info@yafho.com

Wound Dressing


Island Film Dressing

Transparent and waterproof wound dressing with absorbent pad

Features:

1- Transparent 2- Waterproof 3- Breathable 4- Absorbent 5- Hypoallergenic

Ref.	Size(cm)	Package
620206	6x7(3x4)	100pcsx20boxes,2000pcs/ctn
620208	6x10(3x7)	50pcsx24boxes,1200pcs/ctn
620210	9x10(5x6)	30pcsx24boxes,720pcs/ctn
620212	10x12(5x6)	30pcsx24boxes,720pcs/ctn
620215	9x15(5x10)	30pcsx24boxes,720pcs/ctn
620220	10x20(5x15)	25pcsx18boxes,450pcs/ctn
620225	10x25(5x20)	25pcsx18boxes,450pcs/ctn
620230	10x30(5x25)	25pcsx12boxes,300pcs/ctn



Eye Pad(Adhesive)

The Eye bandage is made of non-woven fabric and non-adherent pad. It is specially designed for eye wound care also amblyopia treatment.

Features:

- 1-Adhesive
- 2-Hypoallergenic
- 3-Breathable
- 4-Absorbent
- 5-Other colors available

Ref.	Size(cm)	Package
611302	5.8x8.3(3.2x5)	50pcsx24boxes,1200pcs/ctn
611303	6.5x9.5(3.6x6)	50pcsx24boxes,1200pcs/ctn
612307	5x7.5(2.7x4.5)	100pcsx24boxes,2400pcs/ctn
611307	5x7.5(2.7x4.5)	100pcsx24boxes,2400pcs/ctn
611308	5.5x7.5(3.3x4.9)	100pcsx24boxes,2400pcs/ctn



611302/611303



612307



611307



611308

Yafho
WOUND CARE

Tel:86-20-85295649 Fax:86-20-85295797 <http://www.yafho.com> E-mail:info@yafho.com


I.V. DRESSING

I.V. DRESSING

For fixing PICC/CVC/PIV and small size medical materials.



623178



623115



623113/623114



623252



623273



623277



623419



600328



623251/623271/623275



623272/623270



623509



623309

SECUREFIX +

Transparent Film I.V. Fixation Dressing,
with notch and non woven reinforced.

Ref. (General)	Ref. (HP)	Size(cm)	Package
623252	623252HP	5x5.7	100pcsx20boxes,2000pcs/ctn
623273	623273HP	8.5x10.5	50pcsx20boxes,1000pcs/ctn
623114	623114HP	3.8x4.5	50pcsx24boxes,1200pcs/ctn
623277	623277HP	7x8.5	50pcsx24boxes,1200pcs/ctn
623178	623178HP	8.5x11.5	50pcsx24boxes,1200pcs/ctn
623115	623115HP	6.5x7	50pcsx24boxes,1200pcs/ctn
600328	600328HP	7x9.5	50pcsx24boxes,1200pcs/ctn
623113	623113HP	5x5.7	100pcsx20boxes,2000pcs/ctn
623419	623419HP	7x9	50pcsx24boxes,1200pcs/ctn

HP:W.V.T.R≥1200

General:W.V.T.R≥800

SECUREFIX

Transparent Film I.V. Fixation Dressing

Ref. (General)	Ref. (HP)	Size(cm)	Package
623251	623251HP	4.4x6	100pcsx20boxes,2000pcs/ctn
623279	623279HP	4.4x4.4	100pcsx20boxes,2000pcs/ctn
623271	623271HP	6x7	100pcsx20boxes,2000pcs/ctn
623272	623272HP	6x7	100pcsx20boxes,2000pcs/ctn
623275	623275HP	10x12	50pcsx24boxes,1200pcs/ctn
623270	623270HP	7x9	100pcsx20boxes,2000pcs/ctn
623210	623210HP	9x10	50pcsx24boxes,1200pcs/ctn
623215	623215HP	9x15	50pcsx24boxes,1200pcs/ctn
623220	623220HP	9x20	50pcsx24boxes,1200pcs/ctn
623419	623419HP	9x25	50pcsx24boxes,1200pcs/ctn

HP:W.V.T.R≥1200

General:W.V.T.R≥800

SECUREFIX

I.V. Fixation Dressing ,with pad

Ref.	Size(cm)	Package
623509	6x8(2x2)	100pcsx24boxes,2400pcs/ctn
623309	6x8(2x2)	100pcsx24boxes,2400pcs/ctn

Yafho
WOUND CARE

Tel:86-20-85295649 Fax:86-20-85295797 <http://www.yafho.com> E-mail:info@yafho.com

Advanced Dressing

IV Starter Kit

Our IV Starter Kit was designed to allow clinicians to provide excellent patient care with quality products , while streamlining components ,reducing waste and offering cost effective supply solutions .



Ref.	Package	Composition
623905	50pcsx6boxes,300pcs/ctn	1 I.V Dressing 10x12cm, 1 Catheter Securement Device 9x3.5cm
623906	50pcsx6boxes,300pcs/ctn	2 Alcohol Prep Pad, 1 PVP Prep Pad, 1 I.V Dressing 6.5x7cm, 1 Tourniquet, Latex Free.3/4"x18" 2 Tape 3/4"x18"
623900	50pcsx6boxes,300pcs/ctn	1 I.V Dressing with notch 6x7cm,2 Alcohol Prep Pad, 1 PVP Prep Pad, 4 Surgifix Plaster 2x10cm, 1 Tourbiquet, 1 Dressing Pad 5x5cm, 1 Label Tag.
623902	50pcsx6boxes,300pcs/ctn	1 I.V Dressing 6x7cm, 2 Alcohol Prep Pad, 1 PVP Prep Pad, 4 Surgifix Plaster 2x14cm, 1 Tourbiquet, 1 Dressing Pad 5x5cm, 1 Label Tag.

Wound Visible

Wound Visible's unique and innovative design allows constant monitoring of the wound without the need to disrupt the healing process through unnecessary dressing changes.

Features:

The future of post-operative care is clear to see

Vital protection against serious consequences

Patient Satisfaction

Effective incision site management

Continual monitoring

Help prevent blistering

Ref.	Size(cm)	Package
625208	8x10	20pcsx24boxes,480pcs/ctn
625215	10x15	20pcsx24boxes,480pcs/ctn
625220	10x20	20pcsx24boxes,480pcs/ctn
625225	10x25	20pcsx24boxes,480pcs/ctn
625230	10x30	20pcsx24boxes,480pcs/ctn




Yafho
WOUND CARE

Tel:86-20-85295649 Fax:86-20-85295797 <http://www.yafho.com> E-mail:info@yafho.com

Catheter Securement Device

Ref. (Silicone adhesive)	Ref. (Normal adhesive)	Size(cm)	Package	Notes
600611	600311	16.8x5.7	100pcsx10boxes,1000pcs/ctn	For 2&3 ways Foley catheters sized 12-30 Fr.
600612	600312	16.8x5.7	100pcsx10boxes,1000pcs/ctn	For lines and tubes sized 6-13mm in diameter
600613	600313	9x3.5	50pcsx20boxes,1000pcs/ctn	For arrow PICC/CVC
600614	600314	9x3.5	100pcsx20boxes,2000pcs/ctn	For lines and tubes sized 3-8mm in diameter
600615	600315	7.7x2.2	100pcsx20boxes,2000pcs/ctn	For lines and tubes sized 1.5-4.5mm in diameter
600616	600316	9x3.5	50pcsx20boxes,1000pcs/ctn	For arrow PICC/CVC
600621	600321	4x9.2	30pcsx36boxes,1080pcs/ctn	For lines and tubes
600623	600323	6x7	50pcsx36boxes,1800pcs/ctn	For lines and tubes
600626	600326	6.5x10.5	100pcsx10boxes,1000pcs/ctn	For 2&3 ways Foley catheters sized 12-30 Fr and catheter
600649	600349	9x5(7.5x3.5)	100pcsx10boxes,1000pcs/ctn	For universal femoral and radial catheter
600653	600353	9x5	100pcsx10boxes,1000pcs/ctn	For pediatric foley catheter/tubes sized 12-30Fr in diameter
600610		7x6 5	50pcsx36boxes,1800pcs/ctn	Silicone friendly adhesive, for lines and tubes
	6003371	6 9x8 8	100pcsx8boxes,800pcs /ctn	For lines and tubes
	601106	70x3 5	12pcsx20boxes,240pcs/ctn	Leg band ,non adhesive ,soft material ,supple backing for added skin comfort



Nasal Fix Plaster

Ref.	Size(cm)	Package
600330	6x8	100pcsx20boxes,2000pcs/ctn
600331	8x9	100pcsx20boxes,2000pcs/ctn
600344	6.5x4.5	100pcsx20boxes,2000pcs/ctn
600350	11x4.2	100pcsx10boxes,1000pcs/ctn



Trach- Fix

Tracheostomy Tube Holder

The Tracheostomy Tube Holder is designed to provide secure position and minimize movement of the tracheostomy tube.

Ref.	Size(mm)	Package
601001	430x25	20pcsx20boxes,400pcs/ctn
601002	250x20	30pcsx20boxes,600pcs/ctn
601005	450x20	20pcsx20boxes,400pcs/ctn
601103	450x35	20pcsx20boxes,400pcs/ctn
601102	350x35	20pcsx20boxes,400pcs/ctn



Tracheal Compress

The dressings are specially designed for fixing cannulas. The smooth surface will not stick to the skin. The punched holes in the middle ensure a secure fit of the cannula and tracheal compress.

Ref.	Size(cm)	Package	Notes
687041	10x10	15pcsx30boxes,450pcs/ctn	Foam with split soft
687042	6.5x6.5	15pcsx30boxes,450pcs/ctn	For children made of foam
687094	8.6x7.4	12pcsx30boxes,360pcs/ctn	Foam with split soft
687243	10x10	15pcsx30boxes,450pcs/ctn	Non woven, Double layers with split, soft
687144	10x10	30pcsx30boxes,900pcs/ctn	Non woven with split, soft



Tape

Hypoallergenic Surgical Tape

Sili-Fix

Ref.	Size(cm)	Package	Notes
842250	2.5x500	6rollsx8box,48rolls/ctn	Kind removal Stick repeatedly Skin protectant
842500	5x500	3rollsx8box,24rolls/ctn	

Nonwoven SURGI-FIX

Ref. (Normal)	Ref. (Depend)	Size(cm)	Package
641051	640051	5x1000	100rolls/ctn
641010	640101	10x1000	60rolls/ctn
641151	640151	15x1000	40rolls/ctn
641201	640201	20x1000	40rolls/ctn

TRANSFIX HF





Ref. (Easy Peeling)	Ref.	Size(cm)	Package
6440511	644051	5x1000	100rolls/ctn
6441011	644101	10x1000	60rolls/ctn
6441511	644151	15x1000	40rolls/ctn
6442011	644201	20x1000	40rolls/ctn

Paper Tape,PE Tape

Ref. (Paper Tape)	Ref. (PE Tape)	Size(cm)	Package
642125	643125	1.25x914	24rolls/box,720rolls/ctn
642250	643250	2.5x914	12rolls/box,360rolls/ctn
642500	643500	5.0x914	6rolls/box,180rolls/ctn
642750	643750	7.5x914	6rolls/box,180rolls/ctn

Silk Tape, Zinc Oxide Tape

Ref. (Silk)	Ref. (Zinc)	Size(cm)	Package
645125	647125	1.25x500	24rolls/box,720rolls/ctn
645250	647250	2.5x500	12rolls/box,360rolls/ctn
645500	647500	5x500	6rolls/box,180rolls/ctn
645750	647750	7.5x500	6rolls/box,180rolls/ctn
645100	647100	10x500	6rolls/box,180rolls/ctn



Silicone

Sterile Instrument Isolation Cover

High quality latex free cover provides protection against cross contamination and maintains a sterile field.

Ultrasound Probe Cover(3D Sealing)

Ref.	Size(cm)	Package	Remarks
655283	10x122	60pcsx6boxes,360pcs/ctn	Telescope Folded, Without Gel
655263	10x91	50pcsx6boxes,300pcs/ctn	Telescope Folded, Without Gel
655243	10x61	80pcsx6boxes,480pcs/ctn	Telescope Folded, Without Gel
655343	13x61	80pcsx6boxes,480pcs/ctn	Telescope Folded, Without Gel
655363	13x91	50pcsx6boxes,300pcs/ctn	Telescope Folded, Without Gel
655403	13x122	60pcsx6boxes,360pcs/ctn	Telescope Folded, Without Gel
655284	10x122	30pcsx6boxes,180pcs/ctn	Telescope Folded, With Gel
655264	10x91	30pcsx6boxes,180pcs/ctn	Telescope Folded, With Gel
655244	10x61	30pcsx6boxes,180pcs/ctn	Telescope Folded, With Gel
655344	13x61	30pcsx6boxes,180pcs/ctn	Telescope Folded, With Gel
655364	13x91	30pcsx6boxes,180pcs/ctn	Telescope Folded, With Gel
655404	13x122	30pcsx6boxes,180pcs/ctn	Telescope Folded, With Gel
655566	18x61	80pcsx6boxes,480pcs/ctn	Telescope Folded, Without Gel
655576	18x91	50pcsx6boxes,300pcs/ctn	Telescope Folded, Without Gel
655596	18x122	60pcsx6boxes,360pcs/ctn	Adhesive ,telescope Folded, Without Gel
655567	18x61	30pcsx6boxes,180pcs/ctn	Adhesive ,telescope Folded, With Gel
655577	18x91	30pcsx6boxes,180pcs/ctn	Adhesive ,telescope Folded, With Gel
655597	18x122	30pcsx6boxes,180pcs/ctn	Adhesive ,telescope Folded, With Gel
656000	10x29	100pcsx10boxes,1000pcs/ctn	Ultrasound Probe Condoms Cover

Camera Cover

Ref.	Size(cm)	Package	Remarks
650118	13x242	30pcsx6boxes,180pcs/ctn	telescope folded, without gel
650120	18x244	30pcsx6boxes,180pcs/ctn	telescope folded, without gel
650004	15x250	30pcsx6boxes,180pcs/ctn	telescope folded, without gel

Sterile drapes S3 for SIEMENS Cios Select/Cios Select FD


Part ref	Part number	Ref	Size	Package	usage
650000	Part 1	650910	30cmx210cm	60pcsx2boxes,120pcs /ctn	For C-arm
	Part 2	650512	120cmx60cm	60pcsx2boxes,120pcs /ctn	For Image Intensifier or flat detector
	Part 3	650487	60cmx55cm	60pcsx2boxes,120pcs /ctn	For X-ray tube

Each part can be sold separately

Microscope Cover

Ref.	Size(cm)	Package
650300	120x300	3pcs/box,72pcs/ctn
650250	120x250	3pcs/box,72pcs/ctn

Instrument Cover



Yafho+
WOUND CARE

Tel:86-20-85295649 Fax:86-20-85295797 <http://www.yafho.com> E-mail:info@yafho.com

Tel:86-20-85295649 Fax:86-20-85295797 <http://www.yafho.com> E-mail:info@yafho.com

Yafho+
WOUND CARE


Drape

Incise Drape

The surgical film is made of medical glue and polyurethane film, which is impermeable to bacteria, permeable to oxygen. It is used for protecting the wound site from the bacteria and dirt in an operation.

- Features:
- 1- Waterproof
 - 2- Transparent
 - 3- Elastic
 - 4- Breathable

Ref. (odine)	Ref.	Size(cm)	Package
632455	630455	45x45	20pcs/box,18boxes/ctn
632230	630230	20x30	30pcs/box,18boxes/ctn
632003	630003	28x45	20pcs/box,18boxes/ctn
632005	630005	45x55	10pcs/box,12boxes/ctn
632006	630006	28x30	20pcs/box,12boxes/ctn
632340	630340	30x40	20pcs/box,12boxes/ctn
632098	630098	32x45	20pcs/box,18boxes/ctn
632024	630024	96x60	10pcs/box,10boxes/ctn
632001	630001	10x14	100pcs/box,18boxes/ctn
632182	630182	14x12	100pcs/box,18boxes/ctn
632002	630002	15x28	30pcs/box,18boxes/ctn



Eye Drape

Ref.	Size(cm)	Package
630183	85x65	10pcsx18box,180pcs/ctn

