

I.V. and Catheter Fixation Dressing

--SecureFix

一 次 性 使 用 无 菌 敷 贴

- -留 置 针 固 定 贴

The SecureFix can be used to protect I.V. catheter sites and wounds, to secure devices to the skin and to maintain a moist environment for wound healing.

该敷贴用于保护留置针导管部位和伤口，将设备固定在皮肤上并保持湿润的环境以使伤口愈合

- Provides a waterproof, sterile barrier to external contaminants including liquids, bacteria and viruses.
为伤口提供防水，无菌屏障，阻止外部污染物（包括液体，细菌和病毒）侵入。
- The dressing provides continuous visibility to the I.V. insertion site.
该敷料透明，便于观察留置针使用部位。
- Adhesive is friendly on skin.
粘合剂对皮肤亲和，不易致敏。
- Dressing flexes with skin for greater patient comfort.
敷料可顺应皮肤贴敷，使患者感到舒适。
- High breathable, the skin can breath freely.
敷料透气性好，皮肤可以自由呼吸。

4 Types--四类



SecureFix plus
无纺布加强型



**SecureFix +
catheter
securement
device**
加导管固定型



SecureFix
普通型



SecureFix pad
吸收型

Tel:86-20-85295649 Fax:86-20-85295797 <http://www.yafho.com> E-mail:info@yafho.com

Yafho+
WOUND CARE

单张尺寸：205X280mm

材 料：200克铜版纸

正面

Name/名称	Photo/图片	Ref./型号	Size/尺寸	Application advice/适用建议	Notes/备注说明
SecureFix Plus 无纺布加强型		623178	8.5x11.5cm	CVC/PICC 中心静脉导管 External Jugular Vein 外颈静脉	Border design improve the adhesion, reduce edge-lift 无纺布开窗式设计加强固定 Dot coated border improve the Breathability 无纺布网状涂布, 加强透气 Reduce the incidence of phlebitis caused by i.v. catheters 减少留置导管引起的静脉炎等发生几率
		623115	6.5x7cm	Arterial catheter, PICC, CVC, epidural catheter, PIV (peripheral intravenous catheter) 动脉导管, PICC, CVC, PIV(外周静脉导管)	
		623114	3.8x4.5cm	Peripheral lines - Paediatric 儿童用外周静脉导管	
		623113	5x5.7cm	Peripheral lines - Paediatric 儿童用外周静脉导管	
		623419	7x9cm	Arterial catheter, PICC, CVC, epidural catheter, PIV (peripheral intravenous catheter) 动脉导管, PICC, CVC, PIV(外周静脉导管)	备注Notes: 623419 不含纸框 without paper frame
		623273	8.5x10.5	Peripheral and midline venous catheters, Central venous catheters (including subclavian, jugular, femoral and PICCs), Implanted infusion devices (ports), Dialysis catheters, Pulmonary artery catheters, Arterial catheters 多腔导管、肺动脉导管、血透导管、颈动静脉	Non woven reinforced design improve the adhesion 普通无纺布加强固定 Dot coated border improve the Breathability 无纺布网状涂布, 加强透气 Reduce the incidence of phlebitis caused by i.v. catheters 减少留置导管引起的静脉炎等发生几率
		623277	7x8.5cm	Peripheral lines 周边动脉导管, T-延长管	
		623252	5x5.7cm	Peripheral lines - Paediatric 儿童用外周静脉导管	
		600328	7x9.5cm	For use on peripheral and central venous catheter fixation 周边和中央静脉导管固定	备注Notes: 600328 不含纸框 without paper frame
Securefix +catheter securement device 加导管 固定型		623275 +600313	10X12cm /9x3.5cm	Both secure the i.v. injection site and catheter 适用于留置针及导管的双重固定	The fixation of the catheter reduces the medical risks, such as the retention needle falling off due to accidental pull-out. 导管的固定降低了非预期牵拉导致的留置针脱落等风险
		623271 +600314	6X7cm /9x3.5cm		
Securefix 普通型		623272	6x7cm	Suitable for wing type or with pot retention needle fixation. 适用于翼型或带药壶留置针固定	U notch designed basic type U形凹槽设计基本款 S peeling design on paper frame makes it easy to peel off the backing paper S型纸框撕口设计便于剥离
		623270	7x9cm		
		623251	4.4x6cm	For Pediatric 儿童用	Basic paper frame type 基本纸框款 S-peeling design on paper frame makes it easy to peel off the backing paper S型纸框撕口设计便于剥离
		623271	6x7cm	PIV(peripheral intravenous catheter) 外周静脉注射用PIV	
		623275	10x12cm	CVC/CVP cervical veins, epidural catheters, central veins CVC/CVP 锁骨下静脉, 硬膜外导管, 中心静脉	
Securefix with pad 吸收型		623309	6x8cm	Suitable for wing type or with pot retention needle fixation. 适用于翼型或带药壶留置针固定	U notch designed without frame U形凹槽设计基本款 Absorption pads keep the retention area clean and tidy 吸收垫保持留置部位干净整洁
		623509	6x8cm	Suitable for wing type or with pot retention needle fixation. 适用于翼型或带药壶留置针固定	U notch designed Nonwoven type not waterproof U形凹槽无纺布基本款不防水 Absorption pads keep the retention area clean and tidy 吸收垫保持留置部位干净整洁

Manufacturer: Yafho Bio-Technology Co.,Ltd
生产商：台山市弘毅医疗用品有限公司
Add: No. 12, 9Rd, Changlong Ind. Zone, Taishan city, Guangdong, China
地址：广东江门台山市四九镇长龙工业区九路十二号之一
Tel: 86-020-85295649 Fax: 86-020-85295797
电话：86-20-85295649 传真：86-20-85295797
www.yafho.com

Tel:86-20-85295649 Fax:86-20-85295797 http://www.yafho.com E-mail:info@yafho.com

