

gri[™]

allegeset[®]

DELIVERING PRODUCT SOLUTIONS[™]

ACUTE CARE PRODUCT OFFERING



Product Catalog

ABOUT GRI-ALLESET

DELIVERING PRODUCT SOLUTIONS™

GRI-Alleset is your “one-stop shop” for the widest variety of self-manufactured procedural disposables!

At the same time other companies are rationalizing their product offering, GRI-Alleset approaches each market and its customers with one question in mind... “What is essential to the healthcare provider?”

As a wholly-owned subsidiary of Global Resources International, Inc. (GRI), a market leading medical device manufacturer, the mission of GRI-Alleset is to develop and deliver high-quality, cost-effective product solutions to healthcare providers and supply chain partners as a means to improving the value proposition across the global healthcare marketplace.

Initiating at product ideation and development, through a wide variety of manufacturing capabilities, and offering several cost-effective distribution options, we maintain control of cost and quality throughout the entire product lifecycle.

As an organization we pride ourselves on delivering “true value” - The perfect combination of quality products, competitive pricing, and unparalleled support services.

TABLE OF CONTENTS

Facial Protection

- 2 Eye Shield
- 3 Face Shield

Temperature Maintenance

- 6 Cozee™ Blanket
- 7 Cozee™ Attire
- 9 Temperature
Maintenance

Patient Positioning

- 14 Head Positioners
- 15 Patient Extremities
- 16 Body Positioners
- 17 Patient Positioning
Straps

Prepping

- 20 Sponge Sticks
and Forceps

Surgical Sharps Safety

- 22 Needle Counters
- 27 Magnetic Instrument Pads
- 28 Non-Slip Instrument Pads
- 29 Sharps Transfer Tray and
Scalpel Holder
- 30 Temporary Sharps Holder

Surgical Organization

- 32 Diathermy Instrument
Pouch and Sponge
Counter Bag
- 33 Tube and Cord Holder
and The Finger™

Laparoscopic Accessories

- 36 Veress Needle
- 37 Insufflation Filter and
Tubing Sets
- 39 Antifog Solutions

Surgical Accessories

- 42 Skin Markers
- 44 Electrosurgical Pencils
- 45 Cautery Tip Cleaner
- 46 Towel Clamp
- 47 Light Handle Cover
- 48 Suction Tubing
- 49 Yankauer Suction Tips
- 50 Allesorb Solidifier
- 51 Irrigating Bulb Syringe
& Bulb Syringe
- 52 Cotton Huck Towels
& Low Lint and Non
Woven Towels

Orthopedic Dressings

- 54 Esmark Bandages
- 55 Impervious Stockinette
- 56 Alleflex™ Cohesive
Bandages

Bowls, Basins And Trays

- 58 Basins Liner
- 59 Sponge Basins
- 60 Emesis Basins
- 61 Guidewire Basins
- 62 ECO2 Guidewire Basins
- 63 Lap Tray & Placenta Basins
- 64 Ring Basins
- 65 Rectangular Basins
- 66 Medicine Cups
- 67 Graduated Pitchers
- 68 Assorted Trays

Room Turnover

- 70 Mop Head String
- 71 Mop Head Microfiber
- 72 Room Turnover Kits

Under Pads

- 76 Under Pads

Neonatal And Infant

- 80 Gentleheel® Incision Device
- 81 Baby Blanket
- 82 Baby Cap

Legal Information

FACIAL PROTECTION



EYE SHIELD

The GRI-Alleaset ComfortGuard® Eye Shields are a stylish way to provide complete eye protection during medical and surgical procedures. These protective eye shields offer sleek, updated styling that are comfortable and fit over prescription glasses, thus promoting compliance.

The ComfortGuard® Eye Shields have a special coating to prevent fogging during use, providing for a clear view throughout the procedure.

The lens attachment holes are completely occluded by the frame offering complete protection from splashes and sprays.



- ▶ Fits Over Prescription Glasses - Unique design provides a comfortable and secure fit over prescription glasses so special eye wear is not needed
- ▶ Colored Strips on Lens Top Edge Available in Four Stylish Colors - Offers style and allows for finger-proof handling
- ▶ Reusable Frame - Easily replaceable lenses promote reuse of frame helping to be more cost effective
- ▶ Non-Perforated Lens - Helps prevent fluids from passing through lens attachment holes
- ▶ Innovative Extra Wide Inner Frame Ridge - Fills the gap at the temples, providing better top-down fluid and splash protection
- ▶ Wrap Around Frame - Provides for a more secure fit, preventing slipping, and comfortable placement behind the ears
- ▶ Anti-Fog, Anti-Glare & Anti-Static Clear Lens - Inhibits fog from forming on the lens, provide a glare-free view, and anti-static properties for safety in the operating theater
- ▶ Does not contain natural rubber latex or phthalates - Safe for patients and staff with possible allergies and prevention of leaching of potentially harmful compounds

FACE SHIELD

ACUTE CARE
Facial Protection



GRI-Alleaset continues to expand our line of ComfortGuard® protective apparel products, with not only protection in mind, but also comfort.

The new face shield, a logical extension of our popular ComfortGuard Performance Eye Shields, continues with anti-fog and anti-static coating, for clear, fog-free vision while you perform your most critical work.



In addition to the anti-fog coating, the forehead foam block features slots, so the foam curves around your forehead for the most comfortable fit possible.



The elastic strap is attached with ultrasonic welding. This provides added strap security and eliminates the use of rivets or staples that can easily fall off and get lost.

EYE SHIELD

ACUTE CARE
Facial Protection



Description	REF	
	N Non-Sterile	Case Qty
Eye Shields, Pre-Assembled, Clear with Multicolored Bands	33-0500	200
Eye Shields, Lens, Clear with Multicolored Bands	33-0501	2500
Eye Shields, Frames, Polypropylene, Black	33-0502	200

FACE SHIELD

Description	REF	Case Qty
Face Shield, Full Length, 8.7 in (22 cm) long Slotted Foam, Elastic Band	33-1001	50 / case

TEMPERATURE MAINTENANCE



COZEE™ BLANKETS

The Cozee line of patient warming blankets helps offer superb prevention of hypothermia and maintenance of normothermia through advanced development in fabric technologies and manufacturing practices. Cozee™ Blankets are available in several material and size options to better meet your patient warming needs.



- ▶ Cost Effective Alternative To Active Warming
- ▶ Three-Layer Warmth – 20% Better Insulation Than Cotton Blankets¹
- ▶ High Thermal Retention
- ▶ Non-conductive²
- ▶ Individually Shrink-Wrapped
- ▶ Single Patient Use
- ▶ Latex Free

BLANKETS & WRAPS

Description	Size	REF	
		N Non-Sterile	Case Qty
Warming Blanket, Small	51in x 37in (130cm x 95cm)	68-4115137	40
Warming Blanket, Large	51in x 75in (130cm x 190cm)	68-4115175	20
Reflective Blanket, Small	35in x 40in (88cm x 100cm)	68-480141	60
Reflective Blanket, Medium	48in x 48in (122cm x 122cm)	68-480151	50
Reflective Blanket, Large	48in x 84in (122cm x 213cm)	68-480161	25
Papoose Infant Wrap	-	68-480101	50
Rescue Wrap, Full Enclosure Bag	-	68-480901	10

KITS

Description	Size	REF	
		N Non-Sterile	Case Qty
Kit (Blanket, Adult Cap)	48in x 84in (122cm x 213cm)	68-480961	24
Kit (Blanket, Adult Cap, 2 Leg Warmers)	48in x 48in (122cm x 122cm)	68-480951	20

1. Kansas State University. *A Report on Thermal Insulation Values for Selected Fabrics.* (November 2011). Institute for Environmental Research, Technical Report #12-41.

2. MRI Safety Modeling, LLC. *Report: Evaluation of Conductivity for Cozee Warming Attire Fabric.* (November 2018).

COZEE™ ATTIRE

ACUTE CARE
Temperature Maintenance

Cozee™ helps you meet all your patient care needs by providing you with the following key products in various sizes and configurations—blankets, caps, vests, coats, pants, oximetry mitts, and select kits.

- ▶ High Thermal Retention
- ▶ Allows Easy Access To Surgical Site
- ▶ Cost Effective Alternative To Active Warming
- ▶ Lightweight And Easy To Store
- ▶ Single Patient Use
- ▶ Non-Linting



PATIENT WEAR

Description	REF	
	N Non-Sterile	Case Qty
Cap, Neonatal	68-480221	100
Cap, Infant	68-480231	100
Cap, Child	68-480241	50
Cap, Youth	68-480251	50
Cap, Adult	68-480261	100
Mitt, Oximeter, Infant, Small	68-480441	100
Mitt, Oximeter, Pediatric, Medium	68-480451	100
Mitt, Oximeter, Adult, Large	68-480461	50
Gown, Patient, Silver Front with Green Arms	68-480562	50
Foot Warmers	68-480651	50
Leg Warmers, Adjustable	68-480751	100

CAREGIVER WEAR

Description	REF	
	N Non-Sterile	Case Qty
Jacket, Silver Interior, X-Small	68-480311	50
Jacket, Silver Interior, Small	68-480321	50
Jacket, Silver Interior, Medium	68-480331	50
Jacket, Silver Front with Green Arms, Medium	68-480352	25
Jacket, Silver, Green Interior, Medium	68-480351	25
Jacket, Silver Interior, Large	68-480341	50
Jacket, Silver Front with Green Arms, Large	68-480362	25
Jacket, Silver, Green Interior, Large	68-480361	25
Jacket, Silver, Green Interior, X-Large	68-480371	25
Jacket, Green with Silver Arms, X-Large	68-480372	25
Vest, Silver, Green Interior, Small	68-480841	25
Vest, Silver, Green Interior, Medium	68-480851	25
Vest, Silver, Green Interior, Large	68-480861	25
Vest, Silver, Green Interior, X-Large	68-480871	25
Pants, X-Small	68-480711	50
Pants, Small	68-480721	50
Pants, Medium	68-480731	50
Pants, Large	68-480741	50



TEMPERATURE MAINTENANCE

All of GRI-Alleset's Temperature Monitoring probes are compatible with YSI400 series monitors. The plugs are designed for fluid resistance and incorporate a slip-free, ergonomic grip for comfort, ease of removal, and a firm grip. A longer wire helps keep the bulk of the plug away from the patient to minimize crowding around the face and to ensure it does not create a skin pressure point. They are individually pouched and sterile (EO), made from components free from natural rubber latex and DEHP. They all have an accuracy of $\pm 0.1^{\circ}\text{C}$ at 37°C and $\pm 0.2^{\circ}\text{C}$ at $25\text{-}45^{\circ}\text{C}$.




Description	REF	Case Qty
Interface Cable for Temperature Probes, 5m, Universal	22-IC400UN	10 / case
Interface Cable for Temperature Probes, 5m, HP/Philips	22-IC400HP	10 / case
Interface Cable for Temperature Probes, 5m, Siemens	22-IC400SM	10 / case
Interface Cable for Temperature Probes, 5m, Mindray	22-IC400MR	10 / case

TEMPERATURE MAINTENANCE

GENERAL PURPOSE AND NASAL PROBES

- ▶ Commonly used rectally, nasally, into the esophagus or nasopharynx cavity
- ▶ Smooth sheath with a rounded tip for atraumatic and comfortable insertion
- ▶ Markings every 5cm help to measure insertion depth for safety and patient comfort
- ▶ Available in a soft sheath for less traumatic use in the nasal or nasopharynx cavities

Description	REF	
	S Sterile	Case Qty
Soft Sheath 6 French, Blue (no depth markings)	22-GP406-S	50
Regular Sheath 9 French, White	22-GP409-S	50
Regular Sheath 12 French, White	22-GP412-S	50
Soft Sheath 9 French, Blue	22-NP409-S	50



TEMPERATURE MAINTENANCE

SKIN TEMPERATURE SENSORS

- ▶ Reflective adhesive disk deflects any ambient influences for greater temperature accuracy
- ▶ Adhesive is hypoallergenic to ensure no skin irritation
- ▶ Foam attaches securely to the patient's skin
- ▶ Temperature sensor is embedded in thick foam to reduce pressure points on the patient's skin

Description	REF	
	S Sterile	Case Qty
30in (76cm)	22-SK400-S	50



TEMPERATURE MAINTENANCE


ESOPHAGEAL STETHOSCOPIES

- ▶ Inserted orally into the esophagus on anesthetized patients
- ▶ Specially designed acoustic cuff optimizes sound transmission for best audio quality
- ▶ Beveled tube end minimizes trauma to the patient
- ▶ Acoustic cuff is securely bonded to tube to minimize tissue irritation
- ▶ Two-step earpiece connector incorporates standard luer fittings or other stethoscope tubes for universal connections
- ▶ Tube rigidity is optimized for easy insertion with minimal trauma
- ▶ Available in four sizes

Description	REF	
	S Sterile	Case Qty
9 French	22-ES409-S	50
12 French	22-ES412-S	50
18 French	22-ES418-S	50
24 French	22-ES424-S	50



PATIENT POSITIONING



HEAD POSITIONERS®



Description

Head Cradle, Adult
9" x 8" x 4.25"
(23cm x 20cm x 11cm)

Head Cradle, Pediatric
7" x 6.5" x 3.5"
(18cm x 17cm x 8cm)

N

Non-Sterile

REF

Case Qty

50-150101

36

50-150109

36



Description

Head Cradle with Trach Slit, Adult
9" x 8" x 4.25"
(23cm x 20cm x 11cm)

N

Non-Sterile

REF

Case Qty

50-150102

36



Description

Head Cradle, Stacked with Trach Slit,
Adult 8.5" x 8" x 4"
(22cm x 20cm x 10cm)

N

Non-Sterile

REF

Case Qty

50-150104

36



Description

Head Rest, Donut
9" x 2" (23cm x 5 cm)

N

Non-Sterile

REF

Case Qty

50-150105

36



Description

Head Rest, Donut Multiple
9" x 2" (23cm x 5 cm)

N

Non-Sterile

REF

Case Qty

50-150106

36


PATIENT EXTREMITIES



Description

Ulnar Nerve Protectors
15.5" x 6" x 2"
(39cm x 15cm x 5cm)

		REF
	N Non-Sterile	Case Qty
	50-150201	36 pair



Description

Foot/Heel Protector
7" x 7" x 4" in pair
(18cm x 18cm x 10cm)

		REF
	N Non-Sterile	Case Qty
	50-150204	36 pair



Description

Heel and Elbow Protector
9.25" x 9" x 4"
(24cm x 23cm x 10cm)

		REF
	N Non-Sterile	Case Qty
	50-150205	36 pair



Description

Laminectomy Arm Cradle
24" x 5" x 3"
(61cm x 13cm x 7.5cm)

		REF
	N Non-Sterile	Case Qty
	50-150203	12 pair



Description

Arm Cradle
12.5" x 16.5" x 9"
(32cm x 42cm x 23cm)

		REF
	N Non-Sterile	Case Qty
	50-150202	6

BODY POSITIONERS



Description	REF	
	N Non-Sterile	Case Qty
Foot Drop Support 8.5" x 6.75" x 7.5" (22cm x 17cm x 19cm)	50-150501	12



Description	REF	
	N Non-Sterile	Case Qty
Hip Abduction Pillow 16" x 7.75" x 22" (41cm x 20cm x 53cm)	50-150302	6



Description	REF	
	N Non-Sterile	Case Qty
Body Wedge 24.375" x 10.75" x 8.5" (62cm x 27cm x 22cm)	50-150401	7

PATIENT POSITIONING STRAPS

ACUTE CARE
Patient Positioning

GRI-Alleset's patient positioning straps combines comfort and safety in an affordable package. The portfolio covers a versatile array of disposable, single patient use products for general to procedure specific applications. Developed to be used in acute, alternate site and long-term care, our straps are designed to help you feel secure with your patient positioning needs.



PATIENT POSITIONING STRAPS

KNEE AND BODY STRAPS

- ▶ Constructed with durable hook and loop closures and ring attachments
- ▶ Soft, yet strong material utilizes various fabric technologies to meet end user needs
- ▶ Ideal for use with OR tables, patient beds and gurneys



Hook-and-Loop

	N REF	Case Qty
One Piece Design, grey 2" x 81" (5.1cm x 205.7cm)	50-KB1030	100 each
Two Piece Design, white 3" x 47" (7.6cm x 119.4cm)	50-KB1010	50 each



Foam

	N REF	Case Qty
Two Piece Design, tan 4" x 60" (10.2cm x 152.4cm)	50-KB2030	12 sets



Non-Woven with Metal Buckle

	N REF	Case Qty
Two Piece Design, white 3" x 31" (7.6cm x 78.7cm)	50-KB2010	100 each
Two Piece Design, white 5" x 31" (12.7cm x 78.7cm)	50-KB2011	100 each

ARMBOARD STRAPS


- ▶ Made of soft, comfortable material
- ▶ Secures patient arm to the table



Foam


	N REF	Case Qty
One Piece Design, white 15" x 32" (3.81cm x 81.3cm)	50-AB3010	48 each

PREPPING



SPONGE STICKS

- ▶ Wide variety of sponge sticks
- ▶ Ultrasonic welding used to adhere foam to the stick
- ▶ Available in multiple lengths
- ▶ Absorbent foam available in multiple colors, shapes, and sizes
- ▶ Does not contain natural rubber latex, nano-silver additives or phthalates



Description	REF	
	N Non-Sterile	Case Qty
Stick, Sponge, Blue, 6in (15cm)	25-2008	1000
Stick, Sponge, Orange, 8in (20cm)	25-2009	1000
Stick, Sponge, Rounded Top, Blue, 5.5in (14cm)	25-2010	1000
Stick, Sponge, Angled Handle	25-2011	1000
Stick, Sponge, Rounded Top, Blue, 6.3in (16cm)	25-2012	1000
Stick, Sponge, Vaginal, Rounded Narrow, Orange	25-2020	800


SPONGE FORCEPS



- ▶ Used to hold sponges during prep or fluid collection
- ▶ Available in multiple lengths and styles
- ▶ Does not contain natural rubber latex, nano-silver additives or phthalates

Description	REF	
	N Non-Sterile	Case Qty
Forcep, Sponge, Plastic 5.5in (14cm)	25-0540	500
Forcep, Sponge, Plastic 7.5in (19cm)	25-0550	500
Forcep, Sponge, Plastic 9in (23cm)	25-0562	400
Forcep, White, 5in (13cm)	25-0001	2000
Forcep, All Purpose, 5in (13cm)	25-0106	500

SURGICAL SHARPS SAFETY



SELECT NEEDLE COUNTERS



GRI-Allester needle counters are available in a wide variety of configurations comprised of adhesive strips, foam strips, foam blocks, and magnets, with counts ranging from 10 to 100.

- ▶ Specially selected foam does not shrink during sterilization
- ▶ Available in red or yellow to meet specific market requirements
- ▶ Adhesive tabs secure the counter to desired workspace
- ▶ Easy-to-read numbers allow easy count and placement of blades and needles
- ▶ Clear lid available, providing visualization when closed
- ▶ Does not contain natural rubber latex, nano-silver additives or phthalates

Description	N	REF	S	REF
	Non-Sterile	Case Qty	Sterile	Case Qty
Needle Counter, 10-Count, Magnet, Clear Top	30-01001	100	30-01001-S	100
Needle Counter, 10-Count, Foam Block	30-01010	100	30-01010-S	100
Needle Counter, 10-Count, Foam Block, Magnet	30-01037	100	30-01037-S	100
Needle Counter, 15-Count, Foam Strips, Clear Top	30-01518	96	30-01518-S	100
Needle Counter, 15-Count, Foam Strips, Magnet	30-01585	168	30-01585-S	100
Needle Counter, 15-Count, Magnet, Clear Top	30-01514	168	30-01514-S	100
Needle Counter, 20-Count, Foam Block, Magnet	30-02005	100	30-02005-S	100
Needle Counter, 20-Count, Foam Block, Adhesive	30-02010	100	30-02010-S	100
Needle Counter, 20-Count, Double Magnet	30-02055	168	30-02055-S	100


SELECT NEEDLE COUNTERS

ACUTE CARE
Sharps Safety

Description	N REF		S REF	
	Non-Sterile	Case Qty	Sterile	Case Qty
Needle Counter, 30-Count, Double Magnet	30-03007	100	30-03007-S	100
Needle Counter, 30-Count, Foam Strips	30-03088	96	30-03088-S	100
Needle Counter, 30-Count, Foam Strips, Magnet	30-03085	64	30-03085-S	108
Needle Counter, 30-Count, Magnet, Clear Top	30-03015	64	30-03015-S	108
Needle Counter, 40-Count, Foam Block, Adhesive	30-04027	64	30-04027-S	50
Needle Counter, 40-Count, Foam Block, Magnet	30-04037	64	30-04037-S	50
Needle Counter, 60-Count, Foam Strips	30-06088	64	30-06088-S	108
Needle Counter, 60-Count, Double Magnet	30-06055	64	30-06055-S	108
Needle Counter, 80-Count, Double Foam Block	30-08077	64	30-08077-S	50
Needle Counter, 100-Count, Double Magnet, Foam Block	30-10055	64	30-10055-S	50
Needle Counter, 10-Count, Magnet, Clear Top, Yellow	30-01001Y	100	30-01001Y-S	100
Needle Counter, 10-Count, Foam Block, Yellow	30-01010Y	100	30-01010Y-S	100
Needle Counter, 10-Count, Foam Block, Magnet, Yellow	30-01037Y	100	30-01037Y-S	100
Needle Counter, 15-Count, Foam Strips, Clear Top, Yellow	30-01518Y	96	30-01518Y-S	100
Needle Counter, 15-Count, Foam Strips, Magnet, Yellow	30-01585Y	168	30-01585Y-S	100
Needle Counter, 15-Count, Magnet, Clear Top, Yellow	30-01514Y	168	30-01514Y-S	100
Needle Counter, 20-Count, Foam Block, Magnet, Yellow	30-02005Y	100	30-02005Y-S	100
Needle Counter, 20-Count, Foam Block, Adhesive, Yellow	30-02010Y	100	30-02010Y-S	100
Needle Counter, 20-Count, Double Magnet, Yellow	30-02055Y	168	30-02055Y-S	100
Needle Counter, 30-Count, Double Magnet, Yellow	30-03007Y	100	30-03007Y-S	100
Needle Counter, 30-Count, Foam Strips, Yellow	30-03088Y	96	30-03088Y-S	100
Needle Counter, 30-Count, Foam Strips, Magnet, Yellow	30-03085Y	64	30-03085Y-S	108
Needle Counter, 30-Count, Magnet, Clear Top, Yellow	30-03015Y	64	30-03015Y-S	108
Needle Counter, 40-Count, Foam Block, Adhesive, Yellow	30-04027Y	64	30-04027Y-S	50
Needle Counter, 40-Count, Foam Block, Magnet, Yellow	30-04037Y	64	30-04037Y-S	50
Needle Counter, 60-Count, Foam Strips, Yellow	30-06088Y	64	30-06088Y-S	108
Needle Counter, 60-Count, Double Magnet, Yellow	30-06055Y	64	30-06055Y-S	108
Needle Counter, 80-Count, Double Foam Block, Yellow	30-08077Y	64	30-08077Y-S	50
Needle Counter, 100-Count, Double Magnet, Foam Block, Yellow	30-10055Y	64	30-10055Y-S	50

SECURE NEEDLE COUNTERS

Built-in Blade Remover



Our SECURE Needle Counters maintain all of the features of our SELECT Needle Counters plus the following:

- ▶ Built-in blade remover ensures safety
- ▶ Three (3) hinge design help prevent inadvertent opening and separation
- ▶ Passive locking latch engages automatically and securely
- ▶ Interlocking shells prevent sharps from escaping when closed
- ▶ Does not contain natural rubber latex, nano-silver additives or phthalates

Description	N REF		S REF	
	Non-Sterile	Case Qty	Sterile	Case Qty
Needle Counter, 10-Count, Foam Block, Magnet, Blade Remover	33-01037	100	33-01037-S	100
Needle Counter, 15-Count, Foam Strip w/ Magnet Strip, Adhesive, Blade Remover	33-01580	100	33-01580-S	100
Needle Counter, 15-Count, Foam Strips, Magnet, Blade Remover	33-01585	100	33-01585-S	100
Needle Counter, 20-Count, Foam Block, Magnet, Blade Remover	33-02005	100	33-02005-S	80

SECURE NEEDLE COUNTERS

ACUTE CARE
Sharps Safety

Description	N REF		S REF	
	Non-Sterile	Case Qty	Sterile	Case Qty
Needle Counter, 20-Count, Double Magnet, Blade Remover	33-02007	100	33-02007-S	100
Needle Counter, 20-Count, Foam Block, Adhesive, Blade Remover	33-02010	100	33-02010-S	80
Needle Counter, 20-Count, Foam Strip w/ Magnet Strip, Adhesive, Blade Remover	33-02080	100	33-02080-S	100
Needle Counter, 20-Count, Foam Strip, Adhesive, Blade Remover	33-02081	100	33-02081-S	100
Needle Counter, 30-Count, Double Magnet, Blade Remover	33-03007	100	33-03007-S	100
Needle Counter, 30-Count, Foam Strip w/ Magnet Strip, Adhesive, Blade Remover	33-03080	100	33-03080-S	108
Needle Counter, 30-Count, Foam Strip, Magnet, Blade Remover	33-03085	100	33-03085-S	108
Needle Counter, 30-Count, Foam Strips, Blade Remover	33-03088	100	33-03088-S	100
Needle Counter, 40-Count, Foam Block, Adhesive, Blade Remover	33-04010	100	33-04010-S	50
Needle Counter, 40-Count, Foam Block, Magnet, Blade Remover	33-04037	100	33-04037-S	50
Needle Counter, 60-Count, Double Magnet, Blade Remover	33-06007	100	33-06007-S	108
Needle Counter, 10-Count, Foam Block, Magnet, Blade Remover, Yellow	33-01037Y	100	33-01037Y-S	100
Needle Counter, 15-Count, Foam Strip w/ Magnet Strip, Adhesive, Blade Remover, Yellow	33-01580Y	100	33-01580Y-S	100
Needle Counter, 15-Count, Foam Strips, Magnet, Blade Remover, Yellow	33-01585Y	100	33-01585Y-S	100
Needle Counter, 20-Count, Foam Block, Magnet, Blade Remover, Yellow	33-02005Y	100	33-02005Y-S	80
Needle Counter, 20-Count, Double Magnet, Blade Remover, Yellow	33-02007Y	100	33-02007Y-S	100
Needle Counter, 20-Count, Foam Block, Adhesive, Blade Remover, Yellow	33-02010Y	100	33-02010Y-S	80
Needle Counter, 20-Count, Foam Strip w/ Magnet Strip, Adhesive, Blade Remover, Yellow	33-02080Y	100	33-02080Y-S	100

SECURE NEEDLE COUNTERS

ACUTE CARE
Sharps Safety

Description	N REF		S REF	
	Non-Sterile	Case Qty	Sterile	Case Qty
Needle Counter, 20-Count, Foam Strip, Adhesive, Blade Remover, Yellow	33-02081Y	100	33-02081Y-S	100
Needle Counter, 30-Count, Double Magnet, Blade Remover, Yellow	33-03007Y	100	33-03007Y-S	100
Needle Counter, 30-Count, Foam Strip w/ Magnet Strip, Adhesive, Blade Remover, Yellow	33-03080Y	100	33-03080Y-S	108
Needle Counter, 30-Count, Foam Strip, Magnet, Blade Remover, Yellow	33-03085Y	100	33-03085Y-S	108
Needle Counter, 30-Count, Foam Strips, Blade Remover, Yellow	33-03088Y	100	33-03088Y-S	100
Needle Counter, 40-Count, Foam Block, Adhesive, Blade Remover, Yellow	33-04010Y	100	33-04010Y-S	50
Needle Counter, 40-Count, Foam Block, Magnet, Blade Remover, Yellow	33-04037Y	100	33-04037Y-S	50
Needle Counter, 60-Count, Double Magnet, Blade Remover, Yellow	33-06007Y	100	33-06007Y-S	108

MAGNETIC INSTRUMENT PADS

GRI-Alleset's Magnetic Instrument Pads are designed to safely attract and hold surgical instruments for easy hands-free transferring and safe, quick access to instruments within the sterile field.

- ▶ Encourages safe sharps' passing practices, helping to reduce accidental sharps injuries
- ▶ Powerful raised magnets securely hold instruments in place
- ▶ Non-slip foam backing for enhanced safety
- ▶ Flexible design conforms to shape of patient's body
- ▶ Does not contain natural rubber latex, nano-silver additives or phthalates



Description	N REF		S REF	
	Non-Sterile	Case Qty	Sterile	Case Qty
Instrument Pad, Magnetic, 10x16in (25.4x40.6cm)	31-1016	100	31-1016-S	30
Instrument Pad, Magnetic, 20x16in (50.8x40.6cm)	31-2016	50	31-2016-S	20


NON-SLIP INSTRUMENT PADS

- ▶ Cost-effective solution for hands-free transfer
- ▶ Non-slip surface helps keep instruments in place while also holding securely to the patient drape
- ▶ Encourages safe sharps passing practices, helping to reduce accidental sharps injuries
- ▶ Does not contain natural rubber latex, nano-silver additives or phthalates



Description	N REF		S REF	
	Non-Sterile	Case Qty	Sterile	Case Qty
Instrument Pad, Non-Slip, 20x16in (50.8x40.6cm)	31-11016	100	31-11016-S	108

SHARPS PAD



When compared to the traditional sharps box, GRI-Alleset's Sharps Pad is an economical option to count and dispose of small sharps

- ▶ Resistant to needle and blade rupture
- ▶ Strong adhesive pad
- ▶ Does not contain natural rubber latex, nano-silver additives or phthalates

Description	N REF		S REF	
	Non-Sterile	Case Qty	Sterile	Case Qty
Pad, Sharps	30-2000	360	30-2000-S	100

SHARPS TRANSFER TRAY

ACUTE CARE
Sharps Safety

- ▶ Facilitates the safe transfer of scalpels, syringes, needles and other sharp instruments, eliminating the need for hand-to-hand passing
- ▶ Performs as an Engineering Control in compliance with Bloodborne pathogen regulations such as US OSHA 1910.1030
- ▶ Assists in reducing provider sharps injuries
- ▶ Helps to safely contain and protect delicate instruments when not in use
- ▶ Can be used with a number of sharp instruments, including but not limited to scalpels, syringes, Veress needles, K wires, trocars
- ▶ Not made with natural rubber latex or phthalates



Description	N REF		S REF	
	Non-Sterile	Case Qty	Sterile	Case Qty
Sharps Transfer Tray, Yellow	30-1000	120	30-1000-S	60

SCALPEL HOLDER




- ▶ Accommodates up to 3 scalpels
- ▶ Scalpels maintained in “blade-down” position
- ▶ Adhesive tab allows for secure placement on back table cover or mayo stand
- ▶ Does not contain natural rubber latex, nano-silver additives or phthalates

Description	N REF		S REF	
	Non-Sterile	Case Qty	Sterile	Case Qty
Holder, Scalpel	30-0003	200	30-0003-S	120

TEMPORARY SHARPS HOLDER

GRI-Alleset's Temporary Sharps Holder provides a safe place to hold needles and scalpels before, during, and after surgical procedures.

- ▶ Prevents accidental sharps injuries
- ▶ Adhesive backing allows for secure placement on any work surface
- ▶ Foam block with durable plastic housing
- ▶ Does not contain natural rubber latex, nano-silver additives or phthalates



Description	N REF		S REF	
	Non-Sterile	Case Qty	Sterile	Case Qty
Holder, Sharps, Foam	30-0001	100	30-0001-S	50

SURGICAL ORGANIZATION



DIATHERMY INSTRUMENT POUCH

GRI-Alleset Diathermy Pouches provide secure and sterile storage for instruments while also allowing the instruments to remain in the sterile field during surgery.

- ▶ Three sizes available to meet your needs: 1 pocket, 2 pockets, or 3 pockets
- ▶ Reliable adhesive edge for secure placement
- ▶ Does not contain natural rubber latex, nano-silver additives or phthalates



Description	N REF		S REF	
	Non-Sterile	Case Qty	Sterile	Case Qty
Pouch, Diathermy Instrument, 2 Compartment, 12.6x15.7in (32x40cm)	80-14805	500	80-14805-S	100
Pouch, Diathermy Instrument, 3 Compartment, 11.8x19.7in (30x50cm)	80-14806	500	80-14806-S	100
Pouch, Diathermy Instrument, 1 Compartment, 13x11in (33x28cm)	80-14807	500	80-14807-S	100

SPONGE COUNTER BAG

Clear-backing Sponge Counter Bags help to organize continuous sponge counting and containment in surgery.

- ▶ Easy open pockets for the simple disposal of used sponges
- ▶ Does not contain natural rubber latex, nano-silver additives or phthalates



Description	REF	
	N Non-Sterile	Case Qty
Bag, Sponge Counter	25-2990	250

TUBE AND CORD HOLDER

- ▶ Organizes and securely holds tubes and cords in place
- ▶ Removable hook-and-loop strap allows for added flexibility
- ▶ Adhesive backing fastens easily to drapes
- ▶ Does not contain natural rubber latex, nano-silver additives or phthalates



Description	N REF		S REF	
	Non-Sterile	Case Qty	Sterile	Case Qty
Holder, Tube and Cord, w/ Hook-and-Loop	25-5751	150	25-5751-S	50

THE FINGER™



GRI-Alleset took a fresh-look at the challenge of managing tubing/cables, and created The Finger™. A logical, inexpensive flexible holder that you can place at any position on the surgical drape.


- ▶ One-Piece Design
- ▶ Adjustable Capacity (Up to 3 pieces of suction tubing)
- ▶ Secure (Performs better than 2 piece designs)
- ▶ Easy to Use
- ▶ Cost-Effective (One-piece Design + Automation = Low Cost)
- ▶ Eco-friendly (Less Material = Less Waste)
- ▶ Does not contain natural rubber latex, nano-silver additives or phthalates

Description	REF	
	N Non-Sterile	Case Qty
The Finger™, Tubing and Cable Manager	25-5750	500

LAPAROSCOPIC ACCESSORIES



PERFORMANCE VERESS NEEDLE



The GRI-Alleset Veress Needle, used to establish pneumoperitoneum prior to abdominal endoscopy, has been designed both for patient safety and surgeon control. While the ergonomically designed handle gives the surgeon more feel and control throughout the procedure, the longer beveled edge of the 14-gauge stainless steel needle tip provides for a smoother abdominal entry, thus helping to prevent accidental injury to the patient.

- ▶ Ergonomic Handle – Our unique design provides the surgeon with generous space and raised ridges for improved feel and control
- ▶ Beveled Edge of Needle – Helps facilitate needle entry into the abdominal wall.
- ▶ Red-Safety Indicator – Extending above the plastic handle when the needle tip is exposed, the indicator can be seen clearly even when the surgeon’s fingers cover the handle
- ▶ Blunt-Tip Stylet – Fixed orientation prevents needle tip from blocking airway opening
- ▶ Patented Needle Cap - The patented needle cap is designed to eliminate the creation of particulate during the manufacturing process, by utilizing a snap-on design, as opposed to a slide-on design
- ▶ Not made with natural rubber latex; free of nano-silver additives or phthalates

Description	N REF		S REF	
	Non-Sterile	Case Qty	Sterile	Case Qty
Needle, Veress, 120mm	37-2120	500	37-2120-S	72
Needle, Veress, 150mm	37-2150	500	37-2150-S	72

INSUFFLATION FILTER & TUBING SETS

GRI-Alleset's Insufflation Filter & Tubing Sets prevent bacterial and viral contamination from possible back flow of surgical smoke between patient and insufflator.

- ▶ Bacterial viral filter capable of filtering particles down to 0.1 microns with an average filtration efficiency of 99.986%
- ▶ Luer lock fitting on patient end for connection to Veress needle or primary trocar
- ▶ Flexible, lightweight tubing prevents kinking
- ▶ 10 feet (3 m) long tubing
- ▶ Does not contain natural rubber latex, nano-silver additives or phthalates



Description	N REF		S REF	
	Non-Sterile	Case Qty	Sterile	Case Qty
Tubing, Insufflation w/ Filter, Luer-Lock Fitting	37-1000	100	37-1000-S	25
Tubing, Insufflation w/ Filter, Push Fit Fitting	37-1003	100	37-1003-S	25
Insufflation Tubing Set, 10ft (3m) Rotating Luer Lock Fitting, Color Coded End	37-2000	100	37-2000-S	25
Insufflation Tubing Set, 10ft (3m) Push-On Fitting, Color Coded End	37-2003	100	37-2003-S	25
Insufflation Tubing Set, 10ft (3m) Quick-Connect Fitting, Color Coded End	37-2005	100	37-2005-S	25
Insufflation Tubing 10ft (3m) with 0.1-micron Hi-Flow Filter, 15mm OD Fitting, Color Coded Blue Machine End	37-2015	100	37-2015-S	25
Insufflation Tubing 10ft (3m) with 0.1-micron Hi-Flow Filter, 22mm ID Fitting, Color Coded Blue Machine End	37-2022	100	37-2022-S	25
Insufflation Filter 1ft (30.5cm) Rotating Luer Lock Fitting, Color Coded End	37-2F00	-	37-2F00-S	100
Insufflation Filter, 1ft (30.5cm) Push-On Fitting, Color Coded End	37-2F03	-	37-2F03-S	100
Insufflation Filter, 1ft (30cm) Quick-Connect Fitting, Color Coded End	37-2F05	-	37-2F05-S	100

INSUFFLATION FILTER & TUBING SETS

ACUTE CARE
Lap Accessories

Features	Benefits
Color-coded end	Facilitates aseptic presentation of machine end to circulating nurse
Extra-long (12in / 30cm) machine-end tube	Helps keep filter away from insufflator controls
0.1 Micron filter	Filtering particles down to 0.1 microns with an average efficiency of 99.986% to help prevent infections
Larger filter diameter	Provides up to 71% more surface area for improved air flow that delivers comparable air flow to other companies' "High-flow" tubing
Rotating luer-lock fitting	Patient end of tubing easily connects to Veress needle or primary trocar without twisting tubing
Low-memory tubing	Offers flexibility while helping to prevent kinking
10 foot (3 meter) length	Meets standard market length
Available with five connectors	Machine-end fittings are (1) rotating Luer-lock, (2) push-connect, (3) quick-connect, (4) 15mm OD, or (5) 22mm ID for adaptability to most brands of insufflator machines
Does not contain natural rubber latex, nano-silver additives or phthalates	Safe for patients and staff with possible allergies and prevention of leaching of potentially harmful compounds



2000 Series, Luer Lock



2003 Series, Push-On



2005 Series, Quick-Connect



2015 Series, 15mm OD



2022 Series, 22mm ID

LAPAROSCOPIC ANTI FOG SOLUTION

ACUTE CARE
Lap Accessories

Used by surgeons to reduce fog on laparoscopes, our Anti Fog Solution offers clear results. With a separate 6-gram bottle and sponge in a convenient kit, the sponge can be easily re-saturated as needed throughout the surgical case.

The sponge is a non-abrasive foam that is safe for use on today's delicate scopes and includes a x-ray detectible element for added safety.



Features	Benefits
6-Gram Bottle with Dropper	Plenty of solution to easily moisten the sponge for even the longest surgical procedures
Adhesive Backed Sponge	Allows for convenient and secure placement anywhere on the sterile field
Non-Abrasive Sponge	Reduces possibility of damaging delicate scope lenses
X-Ray Detectible Element in Foam Pad	Assists in x-ray detection if lost in the surgical site
Does not contain natural rubber latex or phthalates	Safe for patients and staff with possible allergies and prevention of leaching of potentially harmful compounds

Description	REF	Case Qty
Anti-Fog Solution 6ml Bottle Non-sterile*	37-001	1100 / case
Sponge Adhesive Pad for Anti-Fog Solution with X-Ray Element Non-sterile	37-002	2000 / case
Anti-Fog Solution Kit 6ml Bottle with Sponge Sterile (gamma)	37-005-S	200 / case

* Non-sterile bottle has been pre-sterilized via gamma irradiation. Liquid inside of bottle is validated as sterile and once assembled and re-packaged, it may be EO sterilized.

SURGICAL ACCESSORIES




SKIN MARKER



GRI-Alleset offers a growing range of Skin Markers for use before and during surgery. Utilizing a special water-based formulation, the Gentian-violet ink is safe for use on the skin. It is non-toxic and will not tattoo like some other markers may do.

The dual-tip marker provides both a fine tip and a standard tip on one convenient skin marker. This dual-tip marker is filled with the same skin-safe Gentian-violet ink.

GRI-Alleset also offers a combo-tip marker, with one end safe for use on skin, while the other “utility” end is formulated with a solvent-based permanent ink for marking plastics, metals, or other items.



- 1 Skin-Safe Ink – Specially formulated ink is non-toxic and will not tattoo
- 2 Available in Standard Tip or Fine Tip - Preferences are met, largely depending on the type of procedure (e.g., plastics versus general use)
- 3 Color-Coded Barrels – Blue for standard tip and white for fine tip makes differentiation simple and clear
- 4 Available Utility Ink – Dual tip marker includes permanent ink for writing on plastic, metal, or other objects
- 5 Does not contain natural rubber latex or phthalates – Safe for patients and staff with possible allergies and prevention of leaching of potentially harmful compounds

Description	N REF		S REF	
	Non-Sterile	Case Qty/ case	Sterile	Case Qty/ case
Combo Tip*	25-D001	1200	25-C001-S	100
Combo Tip* with Ruler	-	-	25-C002-S	100
Combo Tip* with Ruler and Labels	-	-	25-C003-S	100
Standard Tip	-	2000	25-S001-S	100
Standard Tip with Ruler	-	-	25-S002-S	100
Standard Tip with Ruler and Labels	-	-	25-S003-S	100
Fine Tip	-	2000	25-F001-S	100
Fine Tip with Ruler	-	-	25-F002-S	100
Fine Tip with Ruler and Labels	-	-	25-F003-S	100
Dual Tip	-	1800	25-D001-S	100
Dual Tip with Ruler	-	-	25-D002-S	100
Dual Tip with Ruler and Labels	-	-	25-D003-S	100
Standard Tip, "XL" Prep-Resistant	25-SX001	2000	25-SX001-S	100
Standard Tip with Ruler, "XL" Prep-Resistant	-	-	25-SX002-S	100
Standard Tip with Ruler and Labels, "XL" Prep-Resistant	-	-	25-SX003-S	100
Fine Tip, "XL" Prep-Resistant	25-FX001	2000	25-FX001-S	100
Fine Tip with Ruler, "XL" Prep-Resistant	-	-	25-FX002-S	100
Fine Tip with Ruler and Labels, "XL" Prep-Resistant	-	-	25-FX003-S	100

* Combo tip has skin ink in one end and utility ink in the other end.

ELECTROSURGICAL PENCILS

GRI-Alleset offers a standard line of pencils to meet the daily demands of electrosurgery.

GRI-Alleset's electrosurgical pencils are designed with affordability and reliability in mind, with the most common variations on the market readily available. With safety as a primary concern, these pencils have been tested against the most rigorous international electronic safety standard (IEC 60601-1).




- ▶ Economical, yet consistent and reliable functionality
- ▶ Standard three-prong plug fits most of today's electrosurgical units
- ▶ Flexible 3 meter (10 foot) wire has little memory for great drapability
- ▶ Hexagonal collet allows for ideal positional of electrode but without unintended rotation
- ▶ Standard yellow and blue indicators for cut and coagulate functions
- ▶ Switches provide distinct tactical and audible feedback for optimal feel
- ▶ Available with standard or PTFE coated 2cm blades
- ▶ Packaged optionally with holster
- ▶ Not made with natural rubber latex, nano-silver additives or phthalates

Description	N REF		S REF	
	Non-Sterile	Case Qty	Sterile	Case Qty
Push button, 2cm Blade	29-1100	100	29-1100-S	50
Rocker Arm, 2cm Blade	29-1200	100	29-1200-S	50
Push button, 2cm PTFE Coated Blade	29-110C	100	29-110C-S	50
Rocker Arm, 2cm PTFE Coated Blade	29-120C	100	29-120C-S	50
Push button, 2cm Blade, With Holster	-	-	29-1110-S	50
Rocker Arm, 2cm Blade, With Holster	-	-	29-1210-S	50
Push button, 2cm PTFE Coated Blade, With Holster	-	-	29-111C-S	50
Rocker Arm, 2cm PTFE Coated Blade, With Holster	-	-	29-121C-S	50
Holster	29-9000	200	-	-

CAUTERY TIP CLEANERS

ACUTE CARE
Surgical Accessories

- ▶ Effective removal of tissue and debris from cautery tips and blades
- ▶ Coated scratch pad eliminates scraping
- ▶ Adhesive backing for universal placement
- ▶ X-ray detectable
- ▶ Does not contain natural rubber latex, nano-silver additives or phthalates



Description	N REF		S REF	
	Non-Sterile	Case Qty	Sterile	Case Qty
Cleaner, Cautery Tip	25-0200	2000	25-0200-S	100

TOWEL CLAMPS

ACUTE CARE
Surgical Accessories

- ▶ Secure surgical cords, drapes and tubes
- ▶ Easy release tab for repositioning
- ▶ Teeth serrations for added grip
- ▶ Does not contain natural rubber latex, nano-silver additives or phthalates



Description	REF	
	N Non-Sterile	Case Qty
Clamp, Towel, White, 2.5in (6cm)	25-0575	1000

LIGHT HANDLE COVERS


- ▶ Flexible
- ▶ Reduces cost of reprocessing and sterilizing
- ▶ Pushes onto light handle adapter
- ▶ Does not contain natural rubber latex, nano-silver additives or phthalates



Description	N REF		S REF	
	Non-Sterile	Case Qty	Sterile	Case Qty
Bulk	25-3611	1000	-	-
1 Per Pack	-	-	25-3611-S1	120 packs
2 Per Pack	-	-	25-3611-S2	60 packs
3 Per Pack	-	-	25-3611-S3	50 packs

SUCTION TUBING

- ▶ Low coil memory results in greater flexibility
- ▶ Increased wall thickness to help prevent collapsing
- ▶ Color-coded molded female fittings provide easy and secure connections
- ▶ Available in three diameters for various clinical needs:
 - 3/16 in (5mm)
 - 1/4 in (6mm)
 - 9/32 in (7mm)
- ▶ Offered in lengths from 6 to 20 feet (1.8 to 6.1 meters)
- ▶ Presenting two connector options: standard-style and an exclusive EU-style
- ▶ Not made with natural rubber latex and phthalate free for patient safety



Description	N REF		Case Qty	S REF		Case Qty
	Non-Sterile (Standard Connector)	Non-Sterile (EU Connector)		Sterile (Standard Connector)	Sterile (EU Connector)	
3/16 in x 6 ft (5mm x 1.8m)	70-0506C	70-0506EC	50	70-0506C-S	70-0506EC-S	25
3/16 in x 10 ft (5mm x 3.0m)	70-0510C	70-0510EC	50	70-0510C-S	70-0510EC-S	25
3/16 in x 12 ft (5mm x 3.7m)	70-0512C	70-0512EC	50	70-0512C-S	70-0512EC-S	25
3/16 in x 20 ft (5mm x 6.1m)	70-0520C	70-0520EC	50	70-0520C-S	70-0520EC-S	25
1/4 in x 6 ft (6mm x 1.8m)	70-0606C	70-0606EC	50	70-0606C-S	70-0606EC-S	25
1/4 in x 10 ft (6mm x 3.0m)	70-0610C	70-0610EC	50	70-0610C-S	70-0610EC-S	25
1/4 in x 12 ft (6mm x 3.7m)	70-0612C	70-0612EC	50	70-0612C-S	70-0612EC-S	25
1/4 in x 20 ft (6mm x 6.1m)	70-0620C	70-0620EC	50	70-0620C-S	70-0620EC-S	25
9/32 in x 6 ft (7mm x 1.8m)	70-0706C	70-0706EC	50	70-0706C-S	70-0706EC-S	25
9/32 in x 10 ft (7mm x 3.0m)	70-0710C	70-0710EC	50	70-0710C-S	70-0710EC-S	25
9/32 in x 12 ft (7mm x 3.7m)	70-0712C	70-0712EC	50	70-0712C-S	70-0712EC-S	25
9/32 in x 20 ft (7mm x 6.1m)	70-0720C	70-0720EC	50	70-0720C-S	70-0720EC-S	25

YANKAUER SUCTION TIPS

ACUTE CARE
Surgical Accessories



Description (EN)	N REF		S REF	
	Non-Sterile	Case Qty	Sterile	Case Qty
Bulb Tip, Vented	70-50132	130	70-50132-S	50
Bulb Tip, Non-Vented	70-50130	130	70-50130-S	50
Straight Tip, Vented	70-50140	130	70-50140-S	50
Straight Tip, Non-Vented	70-50142	130	70-50142-S	50

ALLESORB SOLIDIFIER

Waste fluid management is an ongoing concern in any healthcare facility. Anything that can help prevent hazardous spills resulting from leaking, cracked, or unsealed suction canisters will be a welcome product.

GRI-Alleaset Allesorb Solidifier is a high-performance super absorbent polymer that, when added to a suction canister, quickly gels the contained liquids to help prevent potentially hazardous spills. Allesorb Solidifier sets the standard in fast and effective waste fluid control.



Features	Benefits
Yorker cap	“Funnel” style cap can be used when only suction ports are available – otherwise, just unscrew the end cap and use the pour spout for faster dispensing
Two Common Sized Bottles	Sized according to the suction canisters you use – no need to measure, just pour in the entire contents for a pre-measured volume
Convenient Dissolvable Pouch	For smaller volumes of waste fluid, just drop-in the pouch. No need to tear and pour: the pouch will dissolve in the liquid
No Shaking Required	Allesorb will gel liquid without the need to shake or stir the waste container, helping to reduce the chances of spills or splashes
Non-Sandy Solid Gel	Allesorb gels in larger clumps, rather than a sandy texture. If a suction canister breaks, the contents of a more solid mass are easier to clean-up
Does not contain natural rubber latex or phthalates	Safe for patients and staff with possible allergies and prevention of leaching of potentially harmful compounds



The dissolvable pouch is a simple and effective. Simply drop the pouch through the suction canister pour spout, and watch as the paper dissolves, simultaneously releasing the solidifier powder and efficiently gelling the liquid.

Description	REF	
	N Non-Sterile	Case Qty
Solidifies 1000cc, Dissolvable Pouch (33gm)	70-S1000	100 / Case
Solidifies 1500cc, Bottle (1.6 oz / 45 gm)	70-S1500	96 / Case
Solidifies 3000cc, Bottle (3.2 oz / 91 gm)	70-S3000	96 / Case

IRRIGATING BULB SYRINGES

- ▶ Pliable bulb for responsive, controllable suction
- ▶ Finger flanges to help prevent roll
- ▶ 60ml
- ▶ Does not contain natural rubber latex, nano-silver additives or phthalates



Description	N REF		S REF	
	Non-Sterile	Case Qty	Sterile	Case Qty
Syringe, Bulb Irrigating, 60ml	25-6700	50	25-6700-S	50

BULB EAR / ULCER SYRINGES

- ▶ Ideal for safely irrigating smaller wounds and secretions
- ▶ Available in 2oz (60ml) or 3oz (89ml)
- ▶ Does not contain natural rubber latex, nano-silver additives or phthalates



Description	N REF		S REF	
	Non-Sterile	Case Qty	Sterile	Case Qty
Syringe, Bulb, 2oz (60ml)	25-6702	200	25-6702-S	20
Syringe, Bulb, 3oz (89ml)	25-6703	150	25-6703-S	20

COTTON HUCK TOWELS

ACUTE CARE
Surgical Accessories

- ▶ Highly absorbent and durable
- ▶ 100% cotton
- ▶ 17x24in (43x61cm)



Description	REF	
	N Non-Sterile	Case Qty
Towel, Cotton, Blue	25-0005	400
Towel, Cotton, White	25-0006	400
Towel, Cotton, Green	25-0007	400
Towel, Cotton, Fenestrated	25-0008	400

LOW LINT NON-WOVEN TOWELS



- ▶ Lower lint alternative to the cotton huck towel
- ▶ Highly absorbent
- ▶ Single-use
- ▶ 17x26in (43x66cm)

Description	REF	
	N Non-Sterile	Case Qty
Towel, Ultra Low Lint, Blue	25-9005	320
Towel, Ultra Low Lint, Green	25-9007	320

ORTHOPEDIC DRESSINGS



ESMARK BANDAGES

GRI-Alleset latex-free bandages offer maximum and consistent pressure to extremities for exsanguination prior to surgery. The smooth finish provides a safe surface to help protect the patient's skin.



- ▶ Advanced washing technology results in minimum powder
- ▶ Does not contain natural rubber latex, non-silver additives or phthalates

Description	N REF		S REF	
	Non-Sterile	Case Qty	Sterile	Case Qty
Bandage, Esmark, Blue, Powderless, 3in x 9ft (7.6 x 275 cm)	25-3090	50	25-3090-S	20
Bandage, Esmark, Blue, Powderless, 3in x 12ft (7.6 x 366 cm)	25-3120	50	25-3120-S	20
Bandage, Esmark, Blue, Powderless, 4in x 9ft (10 x 275 cm)	25-4090	40	25-4090-S	20
Bandage, Esmark, Blue, Powderless, 4in x 12ft (10 x 366 cm)	25-4120	40	25-4120-S	20
Bandage, Esmark, Blue, Powderless, 6in x 9ft (15 x 275 cm)	25-6090	20	25-6090-S	20
Bandage, Esmark, Blue, Powderless, 6in x 12ft (15 x 366 cm)	25-6120	20	25-6120-S	20

IMPERVIOUS STOCKINETTES

GRI-Alleset stockinettes consist of an internal woven material for comfort and an external plastic material that provides an impervious barrier.

- ▶ Available with or without pull tabs and in four sizes to ensure a good fit
- ▶ Does not contain natural rubber latex, nano-silver additives or phthalate



Description	N REF		S REF	
	Non-Sterile	Case Qty	Sterile	Case Qty
Stockinette, Impervious, Small, 7×34in (17.8×86.4cm), without pull tab	85-07340	105	85-07340-S	45
Stockinette, Impervious, Medium, 9×48in (22.9×121.9cm), without pull tab	85-09480	68	85-09480-S	20
Stockinette, Impervious, Large, 12×48in (30.5×121.9cm), without pull tab	85-12480	44	85-12480-S	44
Stockinette, Impervious, X-Large, 16×48in (40.6×121.9cm), without pull tab	85-16480	27	85-16480-S	26
Stockinette, Impervious, Small, 7×34in (17.8×86.4cm), with pull tab	85-07341	84	85-07341-S	30
Stockinette, Impervious, Medium, 9×48in (22.9×121.9cm), with pull tab	85-09481	40	85-09481-S	36
Stockinette, Impervious, Large, 12×48in (30.5×121.9cm), with pull tab	85-12481	44	85-12481-S	40
Stockinette, Impervious, X-Large, 16×48in (40.6×121.9cm), with pull tab	85-16481	27	85-16481-S	26

ALLEGLEX™ COHESIVE BANDAGES

ACUTE CARE
Orthopedic Dressings

- ▶ Cost effective and convenient compared to traditional rubber bandages for post-surgery applications
- ▶ Easy, quick tearing
- ▶ Strong self-adherence with supersoft fabric
- ▶ Industry leading tensile strength
- ▶ Available in 4 convenient sizes
- ▶ Does not contain natural rubber latex, nano-silver additives or phthalates



Description	REF	
	N Non-Sterile	Case Qty
Bandage, Cohesive, Prism, 1in (2.5 cm)	86-10PR	48
Bandage, Cohesive, Prism, 1.5in (3.8 cm)	86-15PR	36
Bandage, Cohesive, Prism, 2in (5 cm)	86-20PR	36
Bandage, Cohesive, Prism, 3in (7.6 cm)	86-30PR	24
Bandage, Cohesive, Prism, 4in (10 cm)	86-40PR	18
Bandage, Cohesive, Kids, 2in (5 cm)	86-20KP	36
Bandage, Cohesive, Tan, 1in (2.5 cm)	86-10TN	48
Bandage, Cohesive, Tan, 1.5in (3.8 cm)	86-15TN	36
Bandage, Cohesive, Tan, 2in (5 cm)	86-20TN	36
Bandage, Cohesive, Tan, 3in (7.6 cm)	86-30TN	24
Bandage, Cohesive, Tan, 4in (10 cm)	86-40TN	18
Bandage, Cohesive, Camo, 1in (2.5 cm)	86-10CA	48
Bandage, Cohesive, Camo, 1.5in (3.8 cm)	86-15CA	36
Bandage, Cohesive, Camo, 2in (5 cm)	86-20CA	36
Bandage, Cohesive, Camo, 3in (7.6 cm)	86-30CA	24
Bandage, Cohesive, Camo, 4in (10 cm)	86-40CA	18
Bandage, Cohesive, Checkered, 1.5in (3.8 cm)	86-15CK	36
Bandage, Cohesive, Checkered, 2in (5 cm)	86-20CK	36
Bandage, Cohesive, Checkered, 3in (7.6 cm)	86-30CK	24

BOWLS, BASINS AND TRAYS



BASIN LINER

GRI-Alleset's Ring Basin Liner is the perfect solution for replacing disposable ring basins, or eliminating the reprocessing of reusable ring basins.

- ▶ Eliminates the acquisition expense and reprocessing of stainless steel ring basins
- ▶ Reduces plastic resin requirement by 20% when compared to disposable ring basins
- ▶ Requires only 19% of the sterile packaging as compared to single, sterile disposable ring basins
- ▶ Requires only 15% of the storage space when compared to single, sterile ring basins



ECO²™ = *ECO*-logical x *ECO*-nomical

Description	Skirt Size	Basin Volume	N REF		S REF	
			Non-Sterile	Case Qty	Sterile	Case Qty
Ring Basin Liner/Basin Stand Cover, Single	52x52in (132x132cm)	6,000cc	80-13506	50	80-13506-S	30
Ring Basin Liner/Basin Stand Cover, Single, Flat	52x52in (132x132 cm)	6,000cc	80-13507	50	-	-
Ring Basin Liner/Basin Stand Cover, Double	69x52in (176x131cm)	6,000cc x 2	80-13509	50	80-13509-S	30

SPONGE BASINS

ACUTE CARE
Basins and Trays

Designed to hold solutions, sponges and other supplies.

- ▶ Graduated in ml, cc and oz for ease of measurement
- ▶ Optional lids available
- ▶ ETO sterilizable
- ▶ Does not contain natural rubber latex, nano-silver additives or phthalates



Description	REF	
	N Non-Sterile	Case Qty
Basin, 150ml, Clear	15-0150	1000
Basin, 250ml	15-0250	600
Basin, 500ml	15-0500	390
Basin, 1000ml	15-1000	340
Basin, 1500ml	15-1500	320
Basin, 2500ml	15-2500	160
Lid Basin, 250ml	15-9025	300
Lid Basin, 1000 & 1500ml	15-9100	160
Lid Basin, 2500ml	15-9250	120

EMESIS BASINS

Multipurpose basin used to hold solutions, instruments and other supplies.

- ▶ Smooth rounded edges
- ▶ Graduated in ml, cc and oz for ease of measurement
- ▶ ETO sterilizable
- ▶ Does not contain natural rubber latex, nano-silver additives or phthalates



700cc

Description	REF	
	N Non-Sterile	Case Qty
Basin, Emesis 700ml	15-0700	320



900cc



Basin, Emesis 800cc
(bowl not included, shown to demonstrate fit only)



Description	REF	
	N Non-Sterile	Case Qty
Basin, Emesis 800ml with Ring Basin arc, Blue	15-0800	TBA
Basin, Emesis 800ml with Ring Basin arc, Clear	15-0800-CL	TBA
Basin, Emesis 800ml with Ring Basin arc, Yellow	15-0800-YL	TBA
Basin, Emesis 900ml, Blue	15-0900	TBA

GUIDEWIRE BASINS

GRI-Alleset’s market-tested designs allow the user allow the user to easily collect the guidewire, while keeping it securely in the basin.

- ▶ Graduated in ml, cc and oz for ease of measurement
- ▶ Optional lid available
- ▶ ETO sterilizable
- ▶ Does not contain natural rubber latex, nano-silver additives or phthalates



Description	REF	
	N Non-Sterile	Case Qty
Basin, Guidewire 2500ml	15-GW2000CT	160
Basin, Guidewire 2500ml, Flat Bottom	15-2000	160
Lid Basin, Guidewire 2500ml	15-9250	120

ECO² GUIDEWIRE BASINS

Introducing the new ECO² Guidewire Bowl.

The ECO² Guidewire Bowl incorporates the basic features as the standard Guidewire Bowls, but with a unique patent pending design that reduces the amount of resin required to produce the bowl, lowers the shipping weight, and improves the nesting in that more bowls can stack in the same height as the legacy design. Consider this image, where 12 of the ECO² Guidewire Bowl stack in the same height as 8 of the legacy bowls.



ECO² = *ECO-logical x ECO-nomical*

- ▶ Less resin required to produce our ECO² Guidewire Bowl: Less plastic = Lower price to you
- ▶ Reduced weight of our ECO² Guidewire Bowl: Reduced freight cost = Lower cost to reach your facility
- ▶ Improved nesting of our ECO² Guidewire Bowl: Same stack height holds 33% more = Reducing your warehouse and wire cart requirements, requires only 15% of the storage space when compared to single, sterile ring basins

Description	REF	
	N Non-Sterile	Case Qty
Basin, Guidewire	16-GW2000CT	240

LAP TRAY

- ▶ Extra-long and narrow design
- ▶ Exclusively designed to contain long laparoscopic and endoscopic instrumentation
- ▶ ETO sterilizable
- ▶ Does not contain natural rubber latex, nano-silver additives or phthalates



Description	REF	
	N Non-Sterile	Case Qty
Laparoscopic Instrument Tray, Long and Narrow, 20 x 56 x 6cm (8 x 22 x 2 in)	15-0822	60

PLACENTA BASINS



Large rectangular basin designed for the collection and inspection of placenta.

- ▶ Graduated in ml, cc and oz for ease of measurement
- ▶ Optional lid available
- ▶ ETO sterilizable
- ▶ Does not contain natural rubber latex, nano-silver additives or phthalates

Description	REF	
	N Non-Sterile	Case Qty
Basin, Placenta 3000ml	15-3000	100
Lid Basin, Placenta 3000ml	15-9300	100

RING BASINS

Designed to hold solutions, instruments and other supplies.

- ▶ Fits all standard ring stands
- ▶ 15-6000 graduated in ml, cc and oz for ease of measurement
- ▶ 15-6500 graduated in cc only
- ▶ ETO sterilizable
- ▶ Does not contain natural rubber latex, nano-silver additives or phthalates



Description	REF	
	N Non-Sterile	Case Qty
Basin Ring, 6000ml	15-6000	60
Basin Ring, 6500ml	15-6500	38

RECTANGULAR BASINS

ACUTE CARE
Basins and Trays



Our largest capacity basin.

- ▶ ETO sterilizable
- ▶ Does not contain natural rubber latex, nano-silver additives or phthalates

Description	REF	
	N Non-Sterile	Case Qty
Basin Rectangular, 8000ml	15-8000	50

MEDICINE CUPS

Highest quality design offers stability in the collection, measurement, and/or distribution of fluids and medications.

- ▶ Various colors available
- ▶ Translucent for easy and accurate measurement
- ▶ Graduated in cc and oz
- ▶ Optional lid available
- ▶ ETO sterilizable
- ▶ Does not contain natural rubber latex, nano-silver additives or phthalates



Description	REF	
	N Non-Sterile	Case Qty
Cup Medicine, 60cc-Clear	15-0060	1150
Cup Medicine, 60cc-Blue	15-0060-BL	1150
Cup Medicine, 60cc-Red	15-0060-RD	1150
Cup Medicine, 60cc-Green	15-0060-GN	1150
Cup Medicine, 60cc-Yellow	15-0060-YL	1150
Lid Cup Medicine, 60cc	15-9006	2000



Description	REF	
	N Non-Sterile	Case Qty
Gallipot, 60cc (2 fl. oz.), clear	15-GP60	750

GRADUATED PITCHERS

Low-profile design for stability and use in procedure packs.

- ▶ Crystal-clear material for easy viewing of contents
- ▶ Pour spout for easy transfer of liquids
- ▶ Integral handle design offers stability
- ▶ Graduated on both sides in cc and oz for easy use by right or left-handed person
- ▶ ETO sterilizable
- ▶ Does not contain natural rubber latex, nano-silver additives or phthalates



Description	REF	
	N Non-Sterile	Case Qty
Pitcher, Graduated 1200cc	15-1200	70
Pitcher, Graduated, 1200cc, ECO2	16-1200	70
Pitcher, Graduated, 1400cc	15-1400	160

 = *ECO-logical x ECO-nomical*

ASSORTED TRAYS



A wide selection of multipurpose trays.

- ▶ Rounded edges prevent tearing or puncture of covers, wraps, and pouches
- ▶ Thick corners prevent cracking or breaking
- ▶ ETO sterilisable
- ▶ Does not contain natural rubber latex, nano-silver additives or phthalates

Description	REF	
	N Non-Sterile	Case Qty
Tray Quarter Shallow 11x13x2.5cm (4.7x5.4x1in)	15-0541	600
Tray Quarter Deep 11x13x5cm (4.7x5.4x2in)	15-0542	600
Tray 2 Compartment 24x11x4.5cm (4.6x9.6x1.8in)	15-0904	300
Tray 3 Compartment 25x24x4.9cm (9.6x10x2in)	15-0910	150
Tray J-Wire 27x22x5cm (9.7x11x2in)	15-0911	150
Tray Half Shallow 24.5x14x5cm (9.6x5.5x2in)	15-0951	300
Tray Half Deep 24.5x14x2.5cm (9.6x5.5x1in)	15-0952	300
Tray Small Platform 28x38x2.5cm (11x15x1in)	15-1115	80
Tray X-Large, X-Deep 45x31.5x6.3cm (17.7x12.5x2.5in)	15-1217	30
Tray Small Mayo 48x31.7x1.5cm (18.9x12.4x0.6in)	15-1219	40
Tray Large Platform 41x38x2.5cm (16.3x14.6x1in)	15-1416	50
Tray Large Mayo 52x38x2cm (16.5x21x0.8in)	15-1622	30

ROOM TURNAROUND



MOP HEAD STRING

GRI-Alleset's range of String Mop Heads, available in cotton, rayon and cotton/rayon blend, are designed with the customer cleaning needs in mind — from healthcare to industrial applications. We understand choosing the right mop head for the job at a price that makes sense is crucial to the buying decision.

- ▶ Available in Cut-End or Looped-End, with or without Fan-Tails
- ▶ Universal headband for use with any style of string mop handles
- ▶ Absorbing over four times their dry weight, they are ideal for tasks like terminal cleaning and more
- ▶ Dries quickly, making them preferred for moisture controlled area applications
- ▶ Extremely low-linting, mitigating lint and debris concerns



Description	REF		Case Qty
	Cut-End	Looped-End	
Cotton 12oz (340g)	MC0112CR12	MC0112LR12	12
Cotton 15oz (425g)	MC0115CR12	MC0115LR12	12
Cotton 15oz with Fan-Tail (425g)	MC0115CF12	MC0115LF12	12
Rayon 12oz (340g)	MR0112CR12	MR0112LR12	12
Rayon 12oz with Fan-Tail (340g)	MR0112CF12	MR0112LF12	12
Rayon 16oz (454g)	MR0116CR12	MR0116LR12	12
Rayon 16oz, with Fan-Tail (454g)	MR0116CF12	MR0116LF12	12
Rayon 18oz (510g)	MR0118CR12	MR0118LR12	12
Rayon 18oz, with Fan-Tail (510g)	MR0118CF12	MR0118LF12	12
Non-Cellulose, 16oz (454g)	50-MOPC01	-	12

MOP HEAD MICROFIBER

When it comes to picking up dirt and dust particles without leaving behind lint, microfiber technology is superior to traditional fibers. This is due to the Micro Looped-Ends, which are created by a split during the manufacturing process.

- ▶ Our mop head surface design decreases drag by approximately 2 times, compared to traditional mops
- ▶ Fitted mounting design slips over the mop handle frame for improved maneuverability and reach
- ▶ This single use product assists in helping reduce cross-contamination by eliminating the need to re-use the product
- ▶ Our microfiber mop cleans effectively, reaching into areas where traditional string mops are too large



Description	REF	
	N Non-Sterile	Case Qty
Mop Head, Flat, Microfiber, 14.6x40cm	MMF0140200	200
Mop Head, Flat, Microfiber, 14.6x45.7cm	MMF0145200	200
Mop Head, Flat, Microfiber, 14.6x51cm	MMF0151200	200

CORE ROOM TURNAROUND KITS®

ACUTE CARE
Room Turnaround

The Right Tools For Faster Turnaround

We've re-imagined how protection, efficiency and speed can be achieved when turning an Operating Room from one procedure to the next. Whatever the procedure, each kit is easily customized to your specific needs. Our Clean Operating Room Environment (CORE) Room Turnaround Kits include components specifically designed to:

- ▶ Minimize cross-contamination risks
- ▶ Increase ease of use
- ▶ Allow your transition team to easily drive down turnaround times



CORE  Room Turnaround Kit®

CUSTOM COMPONENTS

LINENS

- ▶ HydraGuard™ Ultra-Absorbent Impervious Table Covers
- ▶ Patient Lift Sheets
- ▶ Impervious Arm board Covers
- ▶ Impervious Head Rest Covers
- ▶ Impervious Pillow Covers

BAGS

- ▶ Kick Bucket Bags
- ▶ Linen Bags
- ▶ Hazardous Waste Bags
- ▶ Trash Bags


MOP HEAD COVERS

- ▶ Micro Fiber Mop Heads
- ▶ String Mop Heads



CORE ROOM TURNAROUND KITS[®]

ACUTE CARE
Room Turnaround



Super-Absorbent Table Covers

Each CORE Room Turnaround Kit includes a Ultra-Absorbent Table Cover, which utilizes a unique fabric construction designed to wick fluids quickly away from patients, trap them in an absorbent mid-section and prevent strike-through to tables and equipment with an impervious bottom moisture shell that protects patients, staff and equipment.

Stitch-less Sonic Branding

GRI-Alleaset is the first to incorporate sonic bonding into Turnaround Kit linens, maximizing patient, staff, and equipment protection by creating a truly impermeable seam, eliminating holes caused by stitching which can provide a pathway for microbial contamination.

CORE Room Turnaround Kits[®] can be quickly and easily customized specifically for your facility and procedures. Contact one of our specialists to help build your kit.

Turnaround Kit Components may include:

1. Impervious Absorbent Table Cover
2. Lift Sheets
3. Headrest / Pillow Case Covers
4. Armboard Cover w/Strap
5. Microfiber Mop
6. Kick Bucket Bag
7. String Mop
8. Non-Infectious Waste Bag
9. Linen Bag
10. Infectious Waste Bag



CORE ROOM TURNAROUND KITS®

ACUTE CARE
Room Turnaround

CORE Room Turnaround Kits are also available as standard “off-the-shelf” kits.



Standard Turnaround Kits:

Turnaround Kit / REF	49-RTK00	49-RTK01	49-RTK02	49-RTK04
Case Qty	40	36	36	16
Absorbent Table Cover 40X90in (102X229cm)	✓	✓	✓	✓
Lift Sheet 34X66in (86X168cm)	✓	✓	✓	✓
Armboard Covers 13X30in (33X76cm)	✓	—	✓	✓
Headrest Cover 14X15in (36X38cm)	—	—	✓	✓
Drawcord Bag, Clear 38X38in (97X97cm)	—	—	—	✓
Drawcord Bag, Red 38X38in (97X97cm)	—	—	—	✓
Drawcord Bag, Blue 33X38in (84X97cm)	—	—	—	✓
Kickbucket Bag, Red 23X24in (58X61cm)	—	—	—	✓

UNDER PADS



DELIVERING PRODUCT SOLUTIONS™

UNDER PADS

ACUTE CARE
Under Pads

When it comes to Under Pads, it is critical to have a product solution that is not only cost-effective and clinically acceptable, but also addresses the multitude of requirements. With this as our catalyst, GRI-Alleset has developed a comprehensive product solution that sets the industry standards in fluid control applications and is guaranteed to address your patient care needs.

ALLESORB IS AVAILABLE IN THREE VERSIONS:

- ▶ STANDARD allows air-permeability with moderate absorption.
- ▶ SELECT allows air-permeability with high-performance absorption.
- ▶ IMPERVIOUS provides an impermeable layer for maximum protection combined with high performance absorption.



allesorb™
IMPERVIOUS



allesorb®
STANDARD




allesorb™
SELECT


UNDER PADS

ACUTE CARE
Under Pads

Description / Product



Standard



Select



Impervious

TOP LAYER

The top layer quickly wicks fluid from the surface drawing it into the absorbent middle core and away from the patient, staff and equipment.

WICKS

WICKS

WICKS

CORE LAYER

The absorbent middle layer traps and contains the fluids.

ABSORBS
Moderate

ABSORBS
High-Performance

ABSORBS
High-Performance

BASE LAYER

Both the Standard and Select Allesorb Under Pads allow air permeability through the base-layer for the highest standard in patient care, allowing heat to transfer. The base layer helps protect patients, staff and equipment with an impervious moisture barrier.


BREATHS

BREATHS

RETAINS

UNDER PADS

ACUTE CARE
Under Pads

Version	Size	REF	
		 Non-Sterile	Case Qty
Standard	24x36in (61x91cm)	99-UPST2436	60
	31x36in (79x91cm)	99-UPST3136	40
Select	18x24in (46x61cm)	99-UPSE1824	60
	24x36in (61x91cm)	99-UPSE2436	60
	31x36in (79x91cm)	99-UPSE3136	40
Impervious	40x60in (102x150cm)	99-UPIMP4060	30
	40x90in (102x229cm)	99-UPIMP4090	20

NEONATAL AND INFANT




GENTLEHEEL® INCISION DEVICE

ACUTE CARE
Neonatal and Infant

gentleheel® devices offer a safe and efficient way for acquiring blood samples from the heels of newborns and toddlers. Sterile and latex-free, gentleheel® devices are engineered with features that help ensure ease of use and minimal discomfort.

- ▶ Color coded for easier selection
- ▶ Ergonomic design for improved handling
- ▶ Multiple visual indicators reinforce incision site
- ▶ Incision devices are less traumatic than punctures
- ▶ Tri-bevel blade generates a smoother, less traumatic incision
- ▶ Hardened stainless steel blade improves incision accuracy
- ▶ Permanent blade retraction minimizes possible reuse or injury
- ▶ Welded plastic casing prevents accidental exposure to the blade



Description	Incision: Depth / Width	REF	
		 Sterile	Case Qty
gentleheel®, Toddler, Blue	2.00mm / 3.00mm	TGH10X50	50 EA/PK (10 PK/CS)
gentleheel®, Newborn, Green	1.00mm / 2.50mm	NGH10X50	50 EA/PK (10 PK/CS)
gentleheel®, Premie, Lavender	0.85mm / 1.75mm	PGH10X50	50 EA/PK (10 PK/CS)
gentleheel®, Micro-Premie, Yellow	0.65mm / 1.40mm	MPGH10X50	50 EA/PK (10 PK/CS)

BABY BLANKET

ACUTE CARE
Neonatal and Infant

- ▶ Knitted and 100% cotton – disposable but the look & feel of standard blankets
- ▶ Premium soft-napped fabric for enhanced softness and comfort
- ▶ Gender-neutral print options
- ▶ Multiple applications from cribs, general nursing to receiving blanket use
- ▶ Does not contain natural rubber latex, nano-silver additives or phthalates



Description	REF	
	N Non-Sterile	Case Qty
Blanket, Baby, Pink and Blue Stripes, 29.5x40.5in (75x103cm)	25-1340	100
Blanket, Baby, NewB, Pink and Blue Babies, 29.5x40.5in (75x103cm)	25-1341	100
Blanket, Baby, Cozee™, 23.6 X 35.4in (60x90cm)	25-1342	25

BABY CAP

ACUTE CARE
Neonatal and Infant

- ▶ Soft cap rests gently and comfortably on baby's head
- ▶ Knitted in a double-ply unique pattern to prevent heat loss
- ▶ Unisex design with blue and pink stripes
- ▶ One size fits all
- ▶ Does not contain natural rubber latex, nano-silver additives or phthalates



Description	REF	
	N Non-Sterile	Case Qty
Cap, Baby, 6x4in (15x10cm)	25-1336	1400

LEGAL INFORMATION

LEGAL DECLARATIONS


The following are trademarks of Global Resources International, Inc. in the U.S. and other countries. The ® symbol refers to a U.S. registration:

- ▶ Alleflex™
- ▶ Alleset®
- ▶ Allesorb™
- ▶ CORE™
- ▶ ComfortGuard®
- ▶ Cozee™
- ▶ DELIVERING PRODUCT SOLUTIONS™
- ▶ ECO2™
- ▶ Gentleheel®
- ▶ GRI™
- ▶ HyrdaGuard®
- ▶ The Finger™

NOTE: These products may not be available for legal sale in every market. Please check with your sales or account representative for availability of your market of interest.

Copyright © 2019 GRI-Alleset. All rights reserved.

FACILITIES GRI-ALLESET



Corporate Office

GRI-Alleset
4142 Industry Way,
Flowery Branch,
Georgia 30542, USA
Tel: 1.678.866.0550
Email:
info@gri-alleset.com

Europe Office

GRI-Alleset Limited BV
Kingsfordweg 151,
1043 GR
Amsterdam, Netherlands
Tel: 31.20.491.9163
Email:
emea-info@gri-alleset.com

India Office

GRI-Alleset Pvt. Ltd.
331 Arun Chambers,
80, Tardeo Road
Mumbai – 400034
Maharashtra, India
Tel: 91.22.23516311
Email:
india-info@gri-alleset.com

Asia Pacific Office

GRI-Alleset Limited
27/F, World Wide House,
19 Des Voeux Road Central,
Hong Kong
Tel: 852.3975.2898
Email:
apac-info@gri-alleset.com

China Office

GRI-Alleset China
1805 North Honggao Road
Xiuzhou Industry Zone
Jiaxing, Zhejiang,
China 314031
Tel: 86.573.83916501
Email:
china-info@gri-alleset.com

Your Local GRI-Alleset Distributor



DELIVERING PRODUCT SOLUTIONS™

www.gri-alleset.com