



## Innovative Pre-emptive Care

The Devon24D is a portable Deep Vein Thrombosis (DVT) Prevention therapy device. It is fully portable and capable of delivering continuous DVT prevention for up to 24 hours while promoting patient mobility. The Devon24D is created for clinical and home care, requiring no set up and is simple to use with its one-button operation.

The Devon24D helps to enhance proper circulation of blood in the body by mimicking the body's natural muscle contractions, thus increasing venous blood flow to the patient's extremities.



Up to 24  
Hours



Simple  
Operation



Portable

### Specifications

Cycle Time: 60 seconds

Holding Time: 10 seconds

Operation Mode: Up to 24hrs

Charging Time: 3hrs

Weight: 0.32kg / 0.705lbs

Dimensions: 75 x 125 x 33 (mm), 2.95 x 4.7 x 1.3 (in)

AC/DC Adapter: Input 110-240 Volts, AC 50/60Hz, 30VA

MAX. DC Output 5V 2A (2 Plugs)

Battery Type: Lithium Battery 3.7v, 2600mAh (rechargeable)

Default Pressure: 50mmHg,  $\pm 5$ mmHg

