




**SIMPLE  
PREMIUM OA  
KNEE BRACE S2**

Low Profile/Superior Quality/Light Weight

## + 3-POINT LEVERAGE SYSTEM

-Simply provide progressive unloading effectively for OA knee patients



## + INDICATIONS

- Mild to moderate unicompartmental osteoarthritis
- Meniscus procedures/repairs
- Medial or lateral compartment

## + FEATURES

- Double hinge adjustment
- One-step closure system
- Built-in aerospace-grade aluminum frame
- Adjustable range of motion
- Non-slip pads
- Unique design
- Latex free & biological evaluation test approved material (ISO10993-5:2009)






Unique Design



One-Step Closure System



Nonslip Pads



Double Hinge Adjustment

## + ADJUSTMENT PARAMETERS

Varus/valgus angle adjustment range: +/-40°  
 Flexion ROM range: 0°, 30°, 45°, 60°, 75°, 90°  
 Extension ROM range: 0°, 5°, 15°, 20°, 25°

## + MODELS

Left	11007
Right	11011

## + SIZE

Free Size

## + OA SOLUTIONS



OA Knee Sleeve X1




OA Knee Brace S1



Premium OA Knee Brace S2





Dual Upright  
OA Knee Brace M1




**SIMPLE  
PREMIUM OA  
KNEE BRACE S2**

**Simple Medical (Shenzhen) Limited**

 519 Wisdom Valley, No.1010 Bulong Road, Minzhi Street, Longhua District, Shenzhen, China. 518109

 [www.simplemedical.cn](http://www.simplemedical.cn)

 +86-755-86038561